














Bed Rails - BC1180/76



SICC_114A



Bed Rails - BC1180/76

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		IHECRLAS	Assist Rails (Pair)	1	1
2		RAILSAELT-NLA	NLA - Assist Rail Assembly - Left Side	1	1
3		RAILSAERT-NL A	NLA - Assist Rail Assembly - Right Side	1	1
4		RLBRACKAEP- NLA	NLA - Bracket, Mounting, Assist Rail - Left Side	1	1
5		RRBRACKAEP- NLA	NLA - Bracket, Mounting, Assist Rail - Right Side	1	1
6		RLBRACKASTR P-NLA	NLA - Punched Bracket, No Guard - Left Hand	1	1
7		RRBRACKASTR P-NLA	NLA - Punched Bracket, No Guard - Right Hand	1	1
8		IHECBAAS	Set of Assist Bars	1	1
9		RAILABLT-NLA	NLA - Assist Bar - Left Side	1	1
10		RAILABRT-NLA	NLA - Assist Bar - Right Side	1	1
11		KITASSISTR1	Kit, Assist Rail Mounting Hardware	1	1
12		1189259	Kit, Knob - Latch	1	1
			Screw, Shoulder Slotted Head (1/2 x 1/2", 3/8-16)	1	1
			Spring, Compression (3/4 x 1 x 3/64")	1	1
			Knob, LATCH-0002 (3/8-16 x 1-3/8")	1	1
13		1186856	Spring, Compression (3/4 x 1 x 3/64")	1	2
14		1165068	Package, End Plug	10	1