


Bed Rails - BE1180/76 - BE2180/76



Bed Rails - BE1180/76 - BE2180/76

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		IHECRLAS	Assist Rails (Pair)	1	1
2		RAILSAELT-NLA	NLA - Assist Rail Assembly - Left Side	1	1
3		RAILSAERT-NL A	NLA - Assist Rail Assembly - Right Side	1	1
4		IHECBAAS	Set of Assist Bars	1	1
5		RAILABLT-NLA	NLA - Assist Bar - Left Side	1	1
6		RAILABRT-NLA	NLA - Assist Bar - Right Side	1	1
7		KITASSISTR1	Kit, Assist Rail Mounting Hardware	1	1
8		RLBRACKAEP-NLA	NLA - Bracket, Mounting, Assist Rail - Left Side	1	1
9		RRBRACKAEP-NLA	NLA - Bracket, Mounting, Assist Rail - Right Side	1	1
10		RLBRACKASTR P-NLA	NLA - Punched Bracket, No Guard - Left Hand	1	1
11		RRBRACKASTR P-NLA	NLA - Punched Bracket, No Guard - Right Hand	1	1
12		IHECRLFUL-NL A	NLA - Full Rails (Pair)	1	1
12		RAILS0080FE1L -NLA	NLA - Full Rail, Left Side	1	1
12		RAILS0080FE1R -NLA	NLA - Full Rail, Right Side	1	1
13		IHECRLPAR-NL A	NLA - Partial Rails (Pair)	1	1
13		RAILS0080PE1L -NLA	NLA - Partial Rail, Left Side	1	1
13		RAILS0080PE1 R-NLA	NLA - Partial Rail, Right Side	1	1
14		KIT66	Kit, Bracket Rail Bumper	1	1
			Bracket, Rail Bumper	1	4
			Nut, Hex, Flange (5/16-18)	1	4
			Washer, Flat (11/32 x 1-1/4 x 1/8)	1	4
			Bolt, Carriage (5/16-18 x 2-1/4")	1	4
			Anti-Rattle, Rail Spindle	1	4
			Side Rail Arm Holder	1	4
15		SRA0001-NLA	NLA - Short Rail Arm, Right/Left	1	1
16		LRA0001L-NLA	NLA - Long Rail Arm, Left	1	1
17		LRA0001R-NLA	NLA - Long Rail Arm, Right	1	1
18		RLNK0001FASS Y-NLA	NLA - Rail Linkage Assembly, Full	1	1
19		RLNK0001FPNT -NLA	NLA - Rail Linkage, Full Rail	1	1
20		RLNK0001PASS Y-NLA	NLA - Rail Linkage Assembly, Partial	1	1
21		RLNK0001PPNT -NLA	NLA - Rail Linkage, Partial Rail	1	1
22	A	HRL0001A	Head End Side Rail Linkage	1	1
23		SPRML324A-NL A	NLA - Spring, Extension	1	1
24		1189259	Kit, Knob - Latch	1	1

Bed Rails - BE1180/76 - BE2180/76

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
			Screw, Shoulder Slotted Head (1/2 x 1/2", 3/8-16)	1	1
			Spring, Compression (3/4 x 1 x 3/64")	1	1
			Knob, LATCH-0002 (3/8-16 x 1-3/8")	1	1
25		1186856	Spring, Compression (3/4 x 1 x 3/64")	1	2
26		1165068	Package, End Plug	10	1

NOTE: A - Order through SFI Medical Solutions.