


Fabric Clothing Protector

9000XT Only



Fabric Clothing Protector 9000XT Only

| Item Number | Note | Part Number | Description | Quantity Per | |
|--|------|-------------|---|--------------|----------|
| | | | | Package | Assembly |
| 1 | A | 1017447-NLA | NLA - Protector, Fabric Clothing | 1 | 2 |
| 1 | B | 1057388-NLA | NLA - Protector, Fabric Clothing (Pair) | 1 | 2 |
| 2 | C,D | 1054770-NLA | NLA - Kit, Clothing Protector Hardware | 1 | 1 |
| 3 | | 1035976-NLA | NLA - Protector, Fabric Clothing with Hardware (Pair) | 1 | 1 |
| NOTE: A - Before 4/95 B - After 3/95 C - Due to a design change, and Clothing Gaurd hardware needed for chairs built before 2/95 will need to order kit number 1054770 so that the chair can be upgraded to the current design. Previous hardware is NLA. D - Kit includes Cord, Cord Lock, and 2 "D" Rings | | | | | |