


Front Riggings - Footrest Support Assembly GT93H Series



SSM_038A

Front Riggings - Footrest Support Assembly GT93H Series

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,B,C	8881194312	Footrest Support Assembly, Hemi, Color - Right	1	1
1	A,B,C	8881194310	Footrest Support Assembly, Hemi, Color - Left	1	1
1	C	1192525	Footrest Support Assembly, Hemi, Black - Right	1	1
1	C	1192526	Footrest Support Assembly, Hemi, Black - Left	1	1
2		1020719	Anti-Rattle, Black	1	1
3		1006490	Tip, Plastic, Black	1	1
4		1025716	Screw, Phillips Pan Head (#10-32 x 11/16")	1	1
5		1026114	Screw, Phillips Pan Head (#10-32 x 1/2")	1	1
6		1110178	Handle, Release, Black - Right	1	1
6		1110179	Handle, Release, Black - Left	1	1
7		1028846	NLA - Rivnut (#10-32)	1	1
8		1025736	Locknut (1/4-20)	1	1
10		47565X057	Bushing (7/32 x 3/8 x 11/32")	1	1
11		1048085	Housing, Cam	1	1
12		1048086	Lever, Cam, Black	1	1

NOTE: A - Includes Items 2 thru 13. Specify 71P Silver Vein Frame Finish when ordering.
B - NOTE: Footrest Support Tube item #13 is not available, order item #1 as a replacement.
C - Item #9 is not available, order item #1 as a replacement.