

I.V. Rod - Recliner

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1074656	Rod Assembly, I.V., Recliner	1	1
2		04701X041	Rod, I.V., Recliner, Black	1	1
3		1076521	Tube, Support, I.V. Rod, Black	1	1
4		1074652	Holder, I.V. Rod, Black	1	1
5		1027683	Knob, Adjustment, Black	1	1
6		16720X004	Screw, Hex Head (1/4-20 x 1-1/2")	2	2
7		1044762	Bolt, Pivot (5/16-24 x 1-7/8")	1	1
8		1080029	Screw, Phillips Pan Head (6-32 x 3/16")	1	1
9		1008882	Locknut (1/4-20)	1	1
10		16348X009-NLA	NLA - Washer, Steel (3/16 x 9/16 x 1/16")	1	1
11		1010398	Washer, Nylon (1/4 x 1/2 x 1/16")	1	1
12		40157X004	Spacer, Nylon (1/4 x 7/16 x 1/8")	1	1
13		60275X000	Locknut (1/4-20)	1	1
15	B	1095497	Decal, Invacare Logo (1.6)	1	1
		1029377	Instruction Sheet	1	1
NOTE: A - Includes items 2 thru 15 and instruction sheet. B - Not Shown					