


# Armrest, T-Style, Adjustable, Removable Desk Length Armpad (HATD)



# Armrest, T-Style, Adjustable, Removable Desk Length Armpad (HATD)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	SP1639973-NLA	NLA - T-Arm Assembly, Removable, Adjustable Height, Desk Length Armpad - Left	1	1
1	A	SP1639974-NLA	NLA - T-Arm Assembly, Removable, Adjustable Height, Desk Length Armpad - Right	1	1
2		SP1543220	Kit, Desk Length Armpad w/ Hardware	1	1
3		SP1601494	Kit, Hardware, Clothing Guard and Armpad	1	1
4		SP0021-20003	Cap, Cylinder, Black (3/4 Dia. x 1/16")	1	4
5		SP1425601	Clothing Guard, Black, T-Arm (HATD / HATF)	1	1
6	B	SP1606666	Support, Armrest, Plastic	1	1
7		SP1589904	Kit, Locking System, T-Arm (HATD / HATF)	1	1
NOTE: A - Includes items 2-7. B - Item 6 was revised 12/01/17. It is recommended to order (2) ea. of part number SP1606666 for wheelchairs manufactured Before 12/01/17 to match.					