


Armrest, T-Style, Adjustable, Removable Full Length Armpad (HATF)



Armrest, T-Style, Adjustable, Removable Full Length Armpad (HATF)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	SP1639975	T-Arm Assembly, Removable, Adjustable Height, Full Length Armpad, Left	1	1
1	A	SP1639976	T-Arm Assembly, Removable, Adjustable Height, Full Length Armpad, Right	1	1
2		SP1543247	Kit, Full Length Armpad w/ Hardware	1	1
3		SP1606075	Kit, Guard, Clothing Hardware, T-Arm	1	1
4		SP0021-20003	Cap, Cylinder, Black (3/4 Dia. x 1/16")	1	4
5		SP1434725	Clothing Guard, Height Adjustable	1	1
6		SP1435388	Double Nut	1	1
7	B	SP1606666	Support, Armrest, Plastic	1	1
8		SP1589904	Kit, Locking System, T-Arm (HATD / HATF)	1	1
NOTE: A - Includes items 2-8. B - Item 7 was revised 12/01/17. It is recommended to order (2) ea. of part number SP1606666 for wheelchairs manufactured Before 12/01/17 to match.					