


# Headrest Support Assembly, Stealth

## CP175/CP275/CP180/CP280



SHEP\_096A

## Headrest Support Assembly, Stealth CP175/CP275/CP180/CP280

| Item Number | Note    | Part Number | Description   | Quantity Per |          |
|-------------|---------|-------------|---|--------------|----------|
|             |         |             |   | Package      | Assembly |
| 1           | A,J     | 1127408-NLA | NLA - Kit, Headrest 10" TWB Comfort Plus Stealth (CP280)        | 1            | 1        |
| 1           | A,J     | 1127409-NLA | NLA - Kit, Headrest 14" TWB Comfort Plus Stealth (CP180)        | 1            | 1        |
| 1           | A,J     | 1127410-NLA | NLA - Kit, Headrest 10" Multi-Axis Comfort Plus Stealth (CP275) | 1            | 1        |
| 1           | A,J     | 1127411-NLA | NLA - Kit, Headrest 14" Multi-Axis Comfort Plus Stealth (CP175) | 1            | 1        |
| 2           | C       | 1124686     | Kit, Headrest Mounting Hardware                                 | 1            | 1        |
| 3           | C,D     | 1127414     | Kit, Link Extension Stealth                                     | 1            | 1        |
| 4           | B,G,H,I | 0273100     | Kit, Headrest Clamp, Lever Style                                | 1            | 1        |
| 5           | E       | 0237422-NLA | NLA - Clamp, Small Bottom Half                                  | 1            | 1        |
| 6           | E       | 0237522-NLA | NLA - Clamp, Small Top Half                                     | 1            | 1        |
| 7           |         | 1072439-NLA | NLA - Knob, Adjustment, Black (3/4" L Threads)                  | 1            | 1        |
| 8           |         | 0288706-NLA | NLA - Screw, Socket, Flat Head (1/4-20 x 1")                    | 1            | 2        |
| 9           |         | 1025736     | Locknut (1/4-20)  | 1            | 2        |
| 10          | F       | 1160332     | Kit, Headrest Clamp, Knob Style                                 | 1            | 1        |
| 11          |         | 1072439-NLA | NLA - Knob, Adjustment, Black (3/4" L Threads)                  | 1            | 1        |
| 12          | E       | 0237522-NLA | NLA - Clamp, Small Top Half                                     | 1            | 1        |
| 13          |         | 1025736     | Locknut (1/4-20)  | 1            | 2        |
| 14          | E       | 0237422-NLA | NLA - Clamp, Small Bottom Half                                  | 1            | 1        |
| 15          |         | 1009838     | Screw, Socket Button (1/4-20 x 3/4")                            | 1            | 2        |

**NOTE:** A - The mounting hardware will contain a lever and star knob tightening device. Disregard the lever that is provided standard by stealth and use the star knob designed for the Formula TRE systems.

B - Includes items 5-9. Used for ABS Plastic Curved Backs and double spreader bar applications as well.

C - Individual components are not available. The Headrest Mounting Hardware Kit was revised to include a Wire Guard Bracket and is used for Conventional and Contoura Backs. The Headrest Mounting Hardware Kit is also backward compatible in the field.

D - The Stealth Link Extension used on the Stealth Comfort Plus TWB models only

E - Clamp is no longer available individually. See Headrest Clamp Mounting Hardware Kit item 4 or 10 depending on Back Style.

F - Includes items 11-15. Used for Metal Curved Backs. Hardware can also be used with an Invacare/Matrx PCS Back.

G - Item 7 Adjustment Knob was replaced with an Adjustment Lever. The Lever is not available as a separate part. The knob Part # 1072439 is interchangeable with the Lever.

H - Hardware can also be used with an Invacare/Matrx Personal Back, Matrx PB Posture Back and Matrx Elite Backs.

I - For Matrx Elite Backs, item #4 Headrest Clamp Kit attaches to an Adapter Plate that is only included with the Adapter Plate Kit model numbers EHRA-S, EHRA-M or EHRA-L. Refer to the Matrx Rehab Seating and Positioning Pricelist for additional information.

J - Stealth Headrest No Longer Available. See page titled "Headrest Support Assembly - Matrx Elan MEHSP-10, 14" for replacement.