


Propulsion Pro 15 km/h



Propulsion Pro 15 km/h

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1638834	Differential 15 km/h - Pro	--	1
2		1638831	Motor 15 km/h - Pro	--	1
3		1636891	Clutch & Brake 15 km/h (Pro)	--	1
4		1505672	Disclutching Lever 113.2 mm	--	1
5		1571584	Guide	--	1
6		1636874	Speed Sensor, cpl.	--	1
7		1636877	Carbon brush 15 km/h (Pro)	--	1