








Storage



Storage

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1505924	Bracket, Front Basket, cpl.	--	1
2		1637915	Bigger front basket	--	1
3		1637910	Front Box Assembly - black	--	1
4		1637913	Lockabel rear storage box	--	1
5		1505727	Rear Basket, cpl.	--	1
6		1505731	Bag, Seat Rear	--	1
7		1505729	Oxygen Bottle Holder (complete with 3 screws)	--	1