


# Motors and Differentials



SCP\_082A

## Motors and Differentials

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1638833	Differential 12 km/h	--	1
1		1636894	Differential - Ultra	--	1
2		1608266	Motor 15 km/h - Pro	--	1
2		1636875	Motor 10 km/h - Ultra	--	1
3		1636891	Clutch & Brake 15 km/h (Pro)	--	1
3		1636892	Clutch & Brake 10 km/h (Ultra)	--	1
4		1505672	Disclutching Lever 113.2 mm	--	1
5		1571584	Guide	--	1
6		1636877	Carbon brush 15 km/h (Pro)	--	1
6		1636878	Carbon brush 10 km/h (Ultra)	--	1

NOTE: A - Used for Pro model