


# Powerbox and Limit Switch Trigger



SCP\_090A

## Powerbox and Limit Switch Trigger

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1556266	Speed Limit Switch Box Comet Ultra 10 km/h	--	1
1		1556281	Speed Limit Switch Box Comet Pro 15 km/h	--	1
2		1608230	Trigger, limit switch	--	1
3		1553399	Controller, Rhino (120A)	--	1
3		1553400	Controller, Rhino (160A HD)	--	1
	A	1581766	Adjustment Gauge, Speed Reduction Switch	--	1

NOTE: A - Not Shown. For mounting the Speed Reduction Switch, the adjustment Gauge part #1581766 is needed also.