


# Footrest and Hardware (CRFTFR)

When choosing Footrests for Tapered and Non-Tapered Frames, Please use the Footplate Kit Reference Chart when making your selections.



# Footrest and Hardware (CRFTFR)

When choosing Footrests for Tapered and Non-Tapered Frames, Please use the Footplate Kit Reference Chart when making your selections.

Refer to Item #1 Footrest Kit Reference Chart when Selecting for Tapered and Non-Tapered Frames

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	Footrest Kit Reference Chart (CRFTFR)	1	1
1	A	1172746	Kit, Footrest (13" W)	1	1
2			Tube, Footrest (13" W)	1	1
3			Spacer, Footrest	1	2
4			Screw, Socket Head, Cap (M5 x .8 x 30mm)	1	2
5			Nut, Lock (M5 x .8)	1	2
6			Washer, Locking Split (#10)	1	2
1	A	1172747	Kit, Footrest (14" W)	1	1
2			Tube, Footrest (14" W)	1	1
3			Spacer, Footrest	1	2
4			Screw, Socket Head, Cap (M5 x .8 x 30mm)	1	2
5			Nut, Lock (M5 x .8)	1	2
6			Washer, Locking Split (#10)	1	2
1	A	1172748	Kit, Footrest (15" W)	1	1
2			Tube, Footrest (15" W)	1	1
3			Spacer, Footrest	1	2
4			Screw, Socket Head, Cap (M5 x .8 x 30mm)	1	2
5			Nut, Lock (M5 x .8)	1	2
6			Washer, Locking Split (#10)	1	2
1	A	1172749	Kit, Footrest (16" W)	1	1
2			Tube, Footrest (16" W)	1	1
3			Spacer, Footrest	1	2
4			Screw, Socket Head, Cap (M5 x .8 x 30mm)	1	2
5			Nut, Lock (M5 x .8)	1	2
6			Washer, Locking Split (#10)	1	2
1	A	1172750	Kit, Footrest (17" W)	1	1
2			Tube, Footrest (17" W)	1	1
3			Spacer, Footrest	1	2
4			Screw, Socket Head, Cap (M5 x .8 x 30mm)	1	2
5			Nut, Lock (M5 x .8)	1	2
6			Washer, Locking Split (#10)	1	2
1	A	1172751	Kit, Footrest (18" W)	1	1
2			Tube, Footrest (18" W)	1	1
3			Spacer, Footrest	1	2
4			Screw, Socket Head, Cap (M5 x .8 x 30mm)	1	2
5			Nut, Lock (M5 x .8)	1	2
6			Washer, Locking Split (#10)	1	2
1	A	1172752	Kit, Footrest (19" W)	1	1
2			Tube, Footrest (19" W)	1	1
3			Spacer, Footrest	1	2
4			Screw, Socket Head, Cap (M5 x .8 x 30mm)	1	2
5			Nut, Lock (M5 x .8)	1	2
6			Washer, Locking Split (#10)	1	2
1	A	1172753	Kit, Footrest (20" W)	1	1

## Footrest and Hardware (CRFTFR)

When choosing Footrests for Tapered and Non-Tapered Frames, Please use the Footplate Kit Reference Chart when making your selections.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
2			Tube, Footrest (20" W)	1	1
3			Spacer, Footrest	1	2
4			Screw, Socket Head, Cap (M5 x .8 x 30mm)	1	2
5			Nut, Lock (M5 x .8)	1	2
6			Washer, Locking Split (#10)	1	2

NOTE: A - Kit includes items 2 - 6.

### Item #1 Footrest Kit Reference Chart (CRFTFR)

	Non-Tapered Frame 75 and 85 Degree Angle	Tapered Frame 75 Degree Tapered 85 Degree Tapered	Footrest On-Center Dimension (See Illustration) 75 and 85 Degree Angle	Footrest On-Center Dimension (See Illustration) 75, 85 Degree Angle Tapered	Plastic Footrest Cover Non-Tapered	Plastic Footrest Cover Tapered
Chair Width	Part Number	Part Number	Dimension (Inches)	Dimension (Inches)	Part Number	Part Number
13" W	1172746	N/A	11-7/8"	N/A	1179175	N/A
14" W	1172747	N/A	12-7/8"	N/A	1179176	N/A
15" W	1172748	N/A	13-7/8"	N/A	1179177	N/A
16" W	1172749	1172747	14-7/8"	12-7/8"	1179178	1179176
17" W	1172750	1172748	15-7/8"	13-7/8"	1179179	1179177
18" W	1172751	1172749	16-7/8"	14-7/8"	1179180	1179178
19" W	1172752	N/A	17-7/8"	N/A	1179181	N/A
20" W	1172753	N/A	18-7/8"	N/A	1179182	N/A