


Back Canes w/o Push Handles

(LKDWN-Standard) Lock-Down, Folding Height/Angle Adjustable Back



Back Canes w/o Push Handles (LKDWN-Standard) Lock-Down, Folding Height/Angle Adjustable Back

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,B	1170299	NLA - Kit, Back Canes, Adjustable w/o Push Handle Hardware (8 - 11" H)	1	1
2			NLA - Plug Button, Dome, Black (3/4")	1	2
3			NLA - Tube, Insert, Adjustable Back (8 - 11" H)	1	2
4			Spring, Torsion w/ Button	1	2
5			Spacer, Flanged, Black	1	2
1	C,A	1170250	NLA - Kit, Back Canes, Adjustable w/o Push Handle Hardware, (10-14"/16-20" H)	1	1
2			NLA - Plug Button, Dome, Black (3/4")	1	2
3			NLA - Tube, Insert, Adjustable Back (10 - 20" H)	1	2
4			Spring, Torsion w/ Button	1	2
5			Spacer, Flanged, Black	1	2
NOTE: A - Includes items 2 - 5 B - 8-11" High Backs C - 10-14", 12-16", 14-18", and 16-20" High Backs					