


Handrims & Hardware - Natural Fit & Natural Fit LT



6 Point Attaching Wheels (Spoke Wheels)

Handrims & Hardware - Natural Fit & Natural Fit LT

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	Item #1 Natural Fit Handrims - 6 pt. Attaching Wheels	1	1
1		See Chart	Item #1 Natural Fit LT Handrims - 6 pt. Attaching Wheels	1	1
2	A	1145679	Kit, Handrim Attaching Hardware - 6 Point	1	1
			Screw, Phillips Pan Head (#10-32 x 1/2")	1	6
			Locknut (#10-32)	1	6
			NLA - Strap, Securing	1	2
	B	1144492	Instruction Sheet, Natural Fit Handrims	1	1
NOTE: A - Includes Screws, Locknuts, and Securing Straps. B - Not Shown					

Handrims & Hardware - Natural Fit & Natural Fit LT

Item #1 Natural Fit Handrims - 6 Point Attaching Wheels (Spoke Wheels)

Handrims sold as 1 each. Includes all Attaching Hardware (item 2) and Instruction Sheet.

Grip Style Natural Fit	24" Wheels	25" Wheels	26" Wheels
Standard	1145546-NLA	1145547	1145548-NLA
Super Grip	1145550-NLA	1145551	1145552

NOTE: The Natural Fit Handrims have 6 points of attachment to mount to a wheel. The 22" Spoke Wheel (CRF22SP) has four points of attachment and cannot accommodate a Natural Fit Handrim.

Item #1 Natural Fit LT Handrims - 6 Point Attaching Wheels (Spoke Wheels)

Handrims sold as 1 each. Includes all Attaching Hardware (item 2) and Instruction Sheet.

Grip Style Natural Fit - LT	24" Wheels	25" Wheels	26" Wheels
Standard	1163646	1163647-NLA	1163648-NLA
Super Grip	1163649	1163650-NLA	1163651

NOTE: The Natural Fit Handrims have 6 points of attachment to mount to a wheel. The 22" Spoke Wheel (CRF22SP) has four points of attachment and cannot accommodate a Natural Fit Handrim.