


# Accessories - Anti-Tip Assembly



## Accessories - Anti-Tip Assembly

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1172625	Kit, Anti-Tip Assembly w/ Hardware (Pair)	1	1
2			Anti-Tipper	1	2
3			Bracket, Anti-Tip, Black	1	2
4			Spacer (5/16 x 7/16 x 15/32")	1	2
5			Screw, Socket, Cap (M8 x 1.25 x 55mm)	1	2
6			NLA - Locknut (M8 x 1.25)	1	2
7			NLA - Screw, Socket Head (#10-32 x 1-1/4")	1	2

NOTE: A - Includes items 2 - 7.