



# Center Mount Rigging Mounting Hardware (CTMT-FJ/CTMT-F)

Partial view of Seat  
Frame Assembly  
shown for reference



SHEP\_256A

CTMT Front Riggings  
shown for reference



## Center Mount Rigging Mounting Hardware (CTMT-FJ/CTMT-F)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102591	CTMT Front Rigging Mounting Hardware, Ultra Low Maxx	1	1
NOTE: A - CTMT Rigging 12-1/4"-15" STF range does not work for 16" Deep on models ULMT, ULMTR and ULMR.					