


Sound Box Assembly After 7/24/14

Perfecto2VAW



Sound Box Assembly After 7/24/14 Perfecto2VAW

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1131249	Filter, Inlet Sound Reduced	1	1
2		1175254-NLA	NLA - Tech Only - Tubing, PVC (1/4" x 1-29/32")	1	1
3		1106675	Coupling, Inlet Filter	1	1
4		1143435	Resonator Assembly	1	1
5		1175253-NLA	NLA - Tech Only - Restrictor, Inlet Filter w/ Counter Bore	1	1
6		1153693	Foam, Resonator - Top	1	1
7		1153694	Foam, Resonator - Bottom	1	1
8		1128404	Screw, Phillips Hex Washer Head (#8-32 x 3/8")	1	7
9		1143473	Soundbox, Horizontal	1	1
10		2001682	Tie Wrap (8" L)	1	5
11		1185399-NLA	NLA - Tech Only - Spring (10 LG x 0.400 OD)	1	1
12		1185758	Tube, Resonator Inlet Preformed	1	1
13		1136943-NLA	NLA - Tech Only - Fan, Exhaust	1	1
14		1153690-NLA	NLA - Tech Only - Transformer Assembly, 20V - PO2AW	1	1
15		1143491	Rivet, Isolator Mount	1	4
16		1143554	Muffler Assembly	1	1
17		1108810	Screw, Hi-Lo Six Lobe Drive (#10-16 x 5/8")	1	5
18		1147917	Mount, Foam Valve 4-Way	1	1
19		1154030	Foam, Soundbox - Vertical	1	1
20		1143474	Gusset, Soundbox	1	2
21		1186169	Heat Exchanger Assembly	1	1
22		1143582	Clamp, Hose Adjustable (3/8")	1	1
23		1153788	Resonator, Inlet	1	1
24		1153707	Tube, Resonator Inlet Preformed	1	1
25		1153699	Foam, Soundbox - Ceiling	1	1
26		1143445-NLA	NLA - Soundbox, Vertical	1	1
27		60127442	Screw, Slot Hex (10-16)	1	2