




Heel Loops and Toe Loops



**Heel Loop
For
Bi-Lateral Contracture
Footrest**



Heel Loop w/ Ankle Strap



**Toe Loops
(For All Footplates)**

SCM_536A

For Adjustable Angle Footplates - AT5543

Footplate Size	Chair Width	Heel Loop Pair Number	Heel Loop w/ Ankle Strap Single Number (Right)	Heel Loop w/ Ankle Strap Single Number (Left)	Heel Loop w/ Ankle Strap Pair Number	Mounting Hardware Pair Number
Medium (4-1/4" Wide)	12-13" Wide or Tapered 16" Wide	See Adjustable Heel Loop page to order	1035908-NLA	1035909-NLA	1035906	1035931
Large (5-1/4" Wide)	14-15" Wide or Tapered 17-18" Wide	See Adjustable Heel Loop page to order	1035919	1035920	1035918	1035931
X-Large (6-1/4" Wide)	16-17" Wide or Tapered 19-24" Wide	See Adjustable Heel Loop page to order	1035932-NLA	1035933-NLA	1035930	1035931
XX-Large (7-1/4" Wide)	18-24" Wide	See Adjustable Heel Loop page to order	1035932-NLA	1035933-NLA	1035930	1035931

Heel Loops and Toe Loops

For Adjustable Angle Footplates w/ AT5543, Bi-Lateral Contracture Footrest

Footplate Size	Chair Width	Heel Loop Pair Number	Heel Loop w/ Ankle Strap Single Number (Right)	Heel Loop w/ Ankle Strap Single Number (Left)	Heel Loop w/ Ankle Strap Pair Number	Heel Loop Mounting Hardware For Bi-Lateral Contracture Footrest Pair Number	Heel Loop w/ Ankle Strap Mounting Hardware For Bi-Lateral Contracture Footrest Pair Number
Small (3-1/4" Wide)	12-13" Wide or Tapered 16" Wide	1018269	1035908-NLA	1035909-NLA	1035906	1036096	1035931
Medium (4-1/4" Wide)	12-13" Wide or Tapered 16" Wide	1018270	1035908-NLA	1035909-NLA	1035906	1036096	1035931
Large (5-1/4" Wide)	14-15" Wide or Tapered 17-18" Wide	1030958	1035919	1035920	1035918	1036096	1035931
X-Large (6-1/4" Wide)	16-17" Wide or Tapered 19-24" Wide	1069037	1035932-NLA	1035933-NLA	1035930	1036096	1035931
XX-Large (7-1/4" Wide)	18-24" Wide	1069037	1035932-NLA	1035933-NLA	1035930	1036096	1035931

For Aluminum - A1562 Footplates

Footplate Size	Chair Width	Heel Loop Pair Number	Heel Loop w/ Ankle Strap Single Number (Right)	Heel Loop w/ Ankle Strap Single Number (Left)	Heel Loop w/ Ankle Strap Pair Number	Mounting Hardware Pair Number
Small (6-3/4" Wide)	14-15" Wide or Tapered 16-17" Wide	See Adjustable Heel Loop page to order	1035908-NLA	1035909-NLA	1035906	1035931
Medium (7-3/4" Wide)	16-24" Wide or Tapered 16-17" Wide	See Adjustable Heel Loop page to order	1035919	1035920	1035918	1035931

Heel Loops and Toe Loops

For Composite - 1651 Footplates

Footplate Size	Chair Width	Heel Loop Pair Number	Heel Loop w/ Ankle Strap Single Number (Right)	Heel Loop w/ Ankle Strap Single Number (Left)	Heel Loop w/ Ankle Strap Pair Number	Mounting Hardware Pair Number
Small (6" Wide)	13-14" Wide or Tapered 16" Wide	See Adjustable Heel Loop page to order	1035919	1035920	1035918	1035931
Medium (7" Wide)	15-16" Wide or Tapered 17-18" Wide	See Adjustable Heel Loop page to order	1035919	1035920	1035918	1035931
Large (8" Wide)	17-24" Wide or Tapered 19-24" Wide	See Adjustable Heel Loop page to order	1035932-NLA	1035933-NLA	1035930	1035931

For Extra Large - A1564 Footplates

Footplate Size	Chair Width	Heel Loop Pair Number	Heel Loop w/ Ankle Strap Single Number (Right)	Heel Loop w/ Ankle Strap Single Number (Left)	Heel Loop w/ Ankle Strap Pair Number	Mounting Hardware Pair Number
Extra Large (7" Wide)	14-15" Wide	See Adjustable Heel Loop page to order	1035919	1035920	1035918	1035931
Extra Large (8" Wide)	16-24" Wide	See Adjustable Heel Loop page to order	1035919	1035920	1035918	1035931

Toe Loops

Toe Loop Single Number	Toe Loop Pair Number	Mounting Hardware Pair Number
04572X025	1037087	1037089

Heel Loops and Toe Loops

Heel Loop w/ Ankle Strap and Toe Loop Attaching Hardware Kits

Part Number	Part Number
1035931	1037089
For Heel Loop w/ Ankle Strap - Includes: 2) Screw, Phillips Pan Head (1/4-20 x 2-1/4") (2) Spacer (1/4 x 11/32 x 1-3/4") (2) Hex Nut w/ Nylon Cap (1/4-20) (2) Flat Washer (1/4 x 5/8 x 1/16")	For Toe Loops - Includes: (4) Screw, Phillips Pan Head (10-32 x 7/8") (4) Washer (1/4 x 3/4 x 3/32") (4) Locknut (10-32)