


# Caster - 6" Front and Rear - Single Sided Fork



SHEP\_189A  
1641729.iso

## Caster - 6" Front and Rear - Single Sided Fork

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102575	Kit, Head Tube Cap w/ Screw, Front - Right	1	1
1		60102577	Kit, Head Tube Cap w/ Screw, Front - Left	1	1
2	A	60102578	Kit, Front Head Tube Assembly - Right	1	1
2	A	60102576	Kit, Front Head Tube Assembly - Left	1	1
3		40117X006	Bearing (5/8")	1	2
4		1110717	Ring, Tolerance	1	2
5		60104599	Kit, Single Sided Fork Assembly w/ Hardware - Left	1	1
6		60104598	Kit, Single Sided Fork Assembly w/ Hardware - Right	1	1
7		60102605	Kit, Wheel w/ Bearings and Hardware, Black Tire and Black Rim (6" x 2")	1	1
8		60102818	Cap, Head Tube, Black	1	1
NOTE: A - Includes the Front Head Tube, Hardware and items 3 and 4. The Front Head Tube contains two plastic flange bearings (not shown) that are press fitted into place and not sold separately					