


## Accessories - Walker Holder Assembly - CSHW

TDX SP2 Power Wheelchair  
shown for reference only



## Accessories - Walker Holder Assembly - CSWH

| Item Number | Note | Part Number | Description                                                         | Quantity Per |          |
|-------------|------|-------------|---------------------------------------------------------------------|--------------|----------|
|             |      |             |                                                                     | Package      | Assembly |
| 1           | A    | 60108348    | Kit, Walker Holder Assembly, Black - Captain Seat                   | 1            | 1        |
| 2           |      |             | Walker Holder Assembly - Captain Seat                               | 1            | 1        |
| 3           |      |             | Screw, Hex Head (1/4-20 x 1-5/8")                                   | 1            | 4        |
| 4           |      |             | Washer (1/4 x 1/2 x 1/16")                                          | 1            | 4        |
| 5           |      |             | Locknut (1/4-20)                                                    | 1            | 4        |
| 6           |      | 60108347    | Kit, Bracket, Walker Hanger, Black (16, 18, 20, 22" W Captain Seat) | 1            | 2        |
| 6           |      | 60108357    | Kit, Bracket, Walker Hanger, Black (24" W Captain Seat)             | 1            | 2        |

NOTE: A - Includes items 2-5