


Arms - Captain Seat



SHEP_284A

Arms - Captain Seat

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	Item #1 Arm Assembly, Adjustable Height	1	1
2		1142487	Upper Arm Assembly, Desk Length (Right/Left)	1	1
3			Arm Pad, Desk Length, Black Urethane (12-1/4" L x 2-5/8" W)	1	1
4			Screw, Socket Head w/ Patch (5/16-18 x 1-3/4")	1	1
5			Washer, Flat (5/16 x 5/8 x 1/16")	1	2
6			Locknut (5/16-18)	1	1
7			Screw, Hex Head w/ Patch (5/16-18 x 1-1/4")	1	1
8			NLA - Nut, Hex (5/16-18)	1	1
9			Tube, Arm, Upper Adjustable Height, Black	1	1
2	A	1140426	Upper Arm Assembly, Full Length (Right/Left)	1	1
3			Arm Pad, Full Length, Urethane, Black (15" L x 2-5/8" W)	1	1
4			Screw, Socket Head w/ Patch (5/16-18 x 1-3/4")	1	1
5			Washer, Flat (5/16 x 5/8 x 1/16")	1	2
6			Locknut (5/16-18)	1	1
7			Screw, Hex Head w/ Patch (5/16-18 x 1-1/4")	1	1
8			NLA - Nut, Hex (5/16-18)	1	1
9			Tube, Arm, Upper Adjustable Height, Black	1	1
3		1139923	Arm Pad, Desk Length, Black Urethane (12-1/4" L x 2-5/8" W)	1	1
3		1139924	Arm Pad, Full Length, Urethane, Black (15" L x 2-5/8" W)	1	1
10		1114910	Knob, Adjustment, Black (5/16-18 x 3/4")	1	1
11		1114914	Handle, Release, Black (5/16-18 x 1/2")	1	1
12	C	1113646	Tube, Arm, Lower Adjustable Width, Black - Right	1	1
12	C	1113647	Tube, Arm, Lower Adjustable Width, Black - Left	1	1
	D	60108449	Kit, Lower Adjustable Width Arm Tube w/ Hardware, Black (Pair)	1	1
			Knob, Adjustment, Black (5/16-18 x 3/4")	1	2
			Handle, Release, Black (5/16-18 x 1/2")	1	2
			Tube, Arm, Lower Adjustable Width, Black - Right	1	1
			Tube, Arm, Lower Adjustable Width, Black - Left	1	1

NOTE: A - Includes items 2-6
C - Item #12 is not sold separately for wheelchairs with Elevate Option.
D - Part #60108449 used for wheelchairs with Non-Elevate Option only. The Lower Arm Tube horizontal tube measurement is longer than an Elevate Optioned wheelchair.

Item #1 Arm Assembly, Adjustable Height

Wheelchair Option	Description	Right Side Part Number	Left Side Part Number
Non-Elevate	Desk Length Arm	1142494	1142495
		See Note(s) A	See Note(s) A
Non-Elevate	Full Length Arm	1142476	1142477
		See Note(s) A	See Note(s) A
Elevate	Desk Length Arm	60108435	60108436
		See Note(s) A, B	See Note(s) A, B

Arms - Captain Seat

Item #1 Arm Assembly, Adjustable Height

Wheelchair Option	Description	Right Side Part Number	Left Side Part Number
Elevate	Full Length Arm	60108437 See Note(s) B, A	60108438 See Note(s) A, B
NOTE: A - Includes items 2-12 B - Item #12 is not sold separately			