

Rollator Bariatric (Four Wheeled) - Model 66550 Before 11/17/17



Rollator Bariatric (Four Wheeled) - Model 66550 Before 11/17/17

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1175907-NLA	NLA - Kit, Right Brake Handle and Brake Assembly	1	1
1		1175908-NLA	NLA - Kit, Left Brake Handle and Brake Assembly	1	1
2		1176369	Kit, Seat w/ Attaching Hardware	1	1
3		1175909-NLA	NLA - Kit, Rear Wheel Assembly	1	2
4		1176366-NLA	NLA - Kit, Front Wheel Assembly w/ Fork	1	2
5		1176367-NLA	NLA - Kit, Back Rest w/ Attaching Hardware	1	1
6		1176368-NLA	NLA - Bag, Underseat Storage	1	1
	A	1142617-NLA	NLA - Pad, Brake	1	1
	A	1176371	Package, Dust Cap, Black	10	1
	B	1141755-NLA	NLA - Kit, Frame Hardware	1	1
NOTE: A - Not Shown B - Not Shown, Includes (2) Anti-Rattles, (2) Push Buttons, (2) Adjustable Knobs, and (2) Brake Cord Holders.					