


Fendered Side Guards (SIDEGRDS-FDR)



Fendered Side Guards (SIDEGRDS-FDR)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	#1. Fendered Side Guard - Complete	1	1
2	C	SP1601477	Kit, Fendered Side Guard, Medium, w/ L Support - Left	1	1
2	C	SP1601478	Kit, Fendered Side Guard, Medium, w/ L Support - Right	1	1
3	C	SP1601479	Kit, Fendered Side Guard, Large, w/ L Support - Left	1	1
3	C	SP1601480	Kit, Fendered Side Guard, Large, w/ L Support - Right	1	1
4	C	SP1601481	Kit, Fendered Side Guard, Extra Large, w/ L Support - Left	1	1
4	C	SP1601482	Kit, Fendered Side Guard, Extra Large, w/ L Support - Right	1	1
5	C	SP1586515	Fendered Side Guard, Medium - Left	1	1
5	C	SP1586516	Fendered Side Guard, Medium - Right	1	1
6	C	SP1586517	Fendered Side Guard, Large - Left	1	1
6	C	SP1586518	Fendered Side Guard, Large - Right	1	1
7	C	SP1584373	Fendered Side Guard, Extra Large - Left	1	1
7	C	SP1584372	Fendered Side Guard, Extra Large - Right	1	1
8		SP1601485	Kit, Hardware, Fendered Side Guard, L Support	1	1
9	B	SP1651135	Fendered Side Guard Support	1	1

NOTE: A - Includes items 2, 3, 4, 5, 6, 7, 8, and 9.
 B - Item #9 is backward compatible for wheelchairs manufactured Before 05/18/18 also.
 C - Refer to Item #1 chart titled "Fendered Side Guard - Complete" for additional information to help identify the correct Fender Side Guard "SIZE" to order. Need to know both Rear Wheel Size and Rear Seat To Floor (RSTF).

#1. Fendered Side Guard - Complete

Rear Wheel Size	Rear Seat To Floor (RSTF)	Fendered Side Guard Size	Left Side	Right Side
22"	14", 15"	XL = Extra Large	SP1651090	SP1651091
22"	16"	L = Large	SP1651088	SP1651089
22"	17"	M = Medium	SP1651086	SP1651087
22"	18"	N/A	N/A	N/A
24"	14"	N/A	N/A	N/A
24"	15", 16"	XL = Extra Large	SP1651090	SP1651091
24"	17", 18"	L = Large	SP1651088	SP1651089
24"	19"	M = Medium	SP1651086	SP1651087
25"	15"	N/A	N/A	N/A
25"	16", 17", 18"	XL = Extra Large	SP1651090	SP1651091
25"	19"	L = Large	SP1651088	SP1651089
25"	20"	M = Medium	SP1651086	SP1651087