


Fork Support, Fixed Frame



Fork Support, Fixed Frame

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	SP1586271-0005	Kit, Fork Support, w/o Accent, Fixed Frame - Left	1	1
1	A	SP1586285-0005	Kit, Fork Support, w/o Accent, Fixed Frame - Right	1	1
1	A	SP1586271-0641	Kit, Fork Support, Blue Accent, Fixed Frame - Left	1	1
1	A	SP1586285-0641	Kit, Fork Support, Blue Accent, Fixed Frame - Right	1	1
1	A	SP1586271-0642	Kit, Fork Support, Purple Accent, Fixed Frame - Left	1	1
1	A	SP1586285-0642	Kit, Fork Support, Purple Accent, Fixed Frame - Right	1	1
2		SP1589193	Caster Housing Cap, Black	1	1
3		SP1585342-0641	Caster Housing Cap, Blue Accent	1	1
3		SP1585342-0642	Caster Housing Cap, Purple Accent	1	1
4		SP1483872	O-Ring (15,6 x 2,4 NITRILE 70 SH)	1	1
5		SP1586287	Kit, Adjustment Fork Support, Fixed Frame	1	1
6		SP1604465	Kit, Adjustment Part Front Wheel (Pair)	2	1
7		SP1586506	Short Metric Hex Key (5 mm)	1	1

NOTE: A - Includes items 2, 3, and 4.