


Wheels - Rear - Service Parts - Standard and Hub Brake (CS)



SCM_148A

Wheels - Rear - Service Parts - Standard and Hub Brake (CS)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		SP1531456	Tire, Pneumatic Right Run Plus (22 x 1") (25-501)	1	1
1		SP0036-30035	Tire, Pneumatic, Grey (22" x 1-3/8")	1	1
1		SP0036-10014	Tire, Solid, Grey (22" x 1-3/8")	1	1
1		SP1585334	Tire, Pneumatic HP Grey and Black (24" x 1")	1	1
1		SP0036-30001	Tire, Pneumatic, Grey (24" x 1-3/8")	1	1
1		SP0036-10013	Tire, Solid, Grey (24" x 1-3/8")	1	1
2		SP1505299	Kit, Aluminum Handrim w/ Hardware (22")	1	1
2		SP1505920	Kit, Anti-Slip Handrim w/ Hardware (22")	1	1
2		SP1505291	Kit, Aluminum Handrim w/ Hardware (24")	1	1
2		SP1505298	Kit, Anti-Slip Handrim w/ Hardware (24")	1	1
3		SP5750029	Tape, Rim (22")	1	1
3		SP0036-40003	Tape, Rim (24")	1	1
4		SP1061214	Tube, Inner (Presta Valve, 22")	1	1
4		SP0036-20007	Tube, Inner (22" x 1-3/8")	1	1
4		SP1032240	Tube, Inner (24" x 1")	1	1
4		SP0036-20001	Tube, Inner (24" x 1-3/8")	1	1
5		SP1573626	Kit, Handrim Hardware, Spoked Wheel, 4 Point	1	1
6		SP1591280	Kit, Quick Release Axle	1	1
6		SP1591281	Kit, Quick Release Axle for Hub Brake	1	1
7		SPN90003	Kit, Reflector, Spoked Wheel	1	1