





Front Riggings - Pivot/Slide Tube w/ Footplate



Front Riggings - Pivot/Slide Tube w/ Footplate

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,D	1111285	Pivot/Slide Tube w/ Aluminum Footplate Assembly, Black - Right	1	1
1	A,D	1111286	Pivot/Slide Tube w/ Aluminum Footplate Assembly, Black - Left	1	1
1	A	1194329	Pivot Slide Tube Assembly with Composite Footplate, Black - Right	1	1
			Pivot/Slide Tube w/ Plug Button, Footrest, Black Footplate, Composite, Black (6" W x 7-53/64" L)	1	1
			Cap, End, Tube, Black (7/8")	1	1
			Spring, Tension	1	1
			Screw, Phillips Pan Head (#8-32 x 3/8")	1	2
			Plug Button, Black (3/4")	1	1
1	A	1194330	Pivot Slide Tube Assembly with Composite Footplate, Black - Left	1	1
			Pivot/Slide Tube w/ Plug Button, Footrest, Black Footplate, Composite, Black (6" W x 7-53/64" L)	1	1
			Cap, End, Tube, Black (7/8")	1	1
			Spring, Tension	1	1
			Screw, Phillips Pan Head (#8-32 x 3/8")	1	2
			Plug Button, Black (3/4")	1	1
8	A	1110311	Pivot/Slide Tube w/ Composite Footplate Assembly, Black - Right	1	1
8	A	1110312	Pivot/Slide Tube w/ Composite Footplate Assembly, Black - Left	1	1
8	A	1110315	Pivot/Slide Tube w/ Composite Footplate Assembly, Black - Right	1	1
8	A	1110316	Pivot/Slide Tube w/ Composite Footplate Assembly, Black - Left	1	1
9	B	1110344	Pivot/Slide Tube w/ Extra Large Footplate Assembly, Black - Right	1	1
9	B	1110345	Pivot/Slide Tube w/ Extra Large Footplate Assembly, Black - Left	1	1
9	B	1110361	Pivot/Slide Tube w/ Small Aluminum Footplate Assembly, Black - Right	1	1
9	B	1110362	Pivot/Slide Tube w/ Small Aluminum Footplate Assembly, Black - Left	1	1
12	B	1110348	Pivot/Slide Tube w/ Extra Large Footplate Assembly, Black - Right	1	1
12	B	1110349	Pivot/Slide Tube w/ Extra Large Footplate Assembly, Black - Left	1	1
12	B	1110365	Pivot/Slide Tube w/ Small Aluminum Footplate Assembly, Black - Right (14")	1	1
12	B	1110366-NLA	Pivot/Slide Tube w/ Small Aluminum Footplate Assembly, Black - Left (14")	1	1
10	C	1036123	Spring, Leaf	1	1

NOTE: A - Includes items 2-7.
 B - Includes items 2, 3, 6, 7, 10 and 11.
 C - Use one or more leaf springs per assembly as required such that footplate may be moved by hand to the up-right position and remain in place.
 D - NOTE: Part numbers changed from Chrome Pivot/Slide Tube to Black as of 5/16/19