

Pivot Slide Tube with Footplate

Footrest Pivot Slide Tube w/ Aluminum Footplate

Prior Model Number	Current Model Number	Chair Width	Frame side	Part Number
T963S, SH93ES, H903 Chrome	T963A, T943S, SH93ES, TH93S Chrome	14 - 22"	Right	1020232-NLA
T963S, SH93ES, H903 Chrome	T963A, T943S, SH93ES, TH93S Chrome	14 - 22"	Left	1020233-NLA
T963S, H903, SH93ES Black	T963A, SH93ES, T943S, TH93S, Black	14 - 22"	Right	1035200
T963S, H903, SH93ES Black	T963A, SH93ES, T943S, TH93S, Black	14 - 22"	Left	1035201
SH93E, T963A Chrome	SH93E, T963A, T943A, TH93A Chrome	14 - 22"	Right	1017243-NLA
SH93E, T963A Chrome	SH93E, T963A, T943A, TH93A Chrome	14 - 22"	Left	1017244-NLA
SH93E, T963A, SH93E Black	SH93E, T943A, T963A, SH93E, TH93A Black	14 - 22"	Right	1035202
SH93E, T963A, SH93E Black	SH93E, T943A, T963A, SH93E, TH93A Black	14 - 22"	Left	1035203
N/A	T63C Black	14 - 22"	Right	1033059
N/A	T63C Black	14 - 22"	Left	1033060
N/A	T63C Chrome	14 - 22"	Right	1017224-NLA
N/A	T63C Chrome	14 - 22"	Left	1017225-NLA
H93E, H93M	H93E, H93M Chrome	14 - 15"	Right	1020232-NLA
H93E, H93M	H93E, H93M Chrome	14 - 15"	Left	1020233-NLA
H93E, H93M	H93E, H93M Chrome	16 - 22"	Right	1028905-NLA
H93E, H93M	H93E, H93M Chrome	16 - 22"	Left	1028906-NLA
H93M w/ FPP	H93M w/ FPP	14 - 15"	Right	1031936-NLA
H93M w/ FPP	H93M w/ FPP	14 - 15"	Left	1031937-NLA
H93M w/ FPP	H93M w/ FPP Chrome	16 - 22"	Right	1033381-NLA
H93M w/ FPP	H93M w/ FPP Chrome	16 - 22"	Left	1033382-NLA
N/A	T93HA Black	18 - 22"	Right	1110344
N/A	T93HA Black	18 - 22"	Left	1110345
N/A	T93HSA, T93HA Black	14 - 15"	Right	1110361

Pivot Slide Tube with Footplate

Footrest Pivot Slide Tube w/ Aluminum Footplate

Prior Model Number	Current Model Number	Chair Width	Frame side	Part Number
N/A	T93HSA, T93HA Black	14 - 15"	Left	1110362
N/A	T93HSA, T93HA Chrome	14 - 15"	Right	1110363-NLA
N/A	T93HSA, T93HA Chrome	14 - 15"	Left	1110364-NLA
N/A	GT93HA/T93HA, T93HAB - Black	16 - 22"	Right	1111285
N/A	GT93HA/T93HA, T93HAB - Black	16 - 22"	Left	1111286
N/A	T93HA, T93HAB - Black	16 - 22"	Right	1111287-NLA
N/A	T93HA, T93HAB - Black	16 - 22"	Left	1111288-NLA
NOTE: 2" and 4" Longer Pivot Slide Tube with Footplate Assembly is not available.				

Pivot Slide Tube with Footplate

Legrest Pivot Slide Tube w/ Aluminum Footplate

Prior Model #	Current Model #	Chair Width	Side	Part Number
SH904MS	SH904MS, T94HS Chrome	14 - 15"	Right	1020234-NLA
SH904MS	SH904MS, T94HS Chrome	14 - 15"	Left	1020235-NLA
SH904M	SH904M	16 - 18"	Right	1017237-NLA
SH904M	SH904M	16 - 18"	Left	1017238-NLA
H904	T944S	14 - 15"	Right	1020234-NLA
H904	T944S	14 - 15"	Left	1020235-NLA
N/A	T944A, TH94A Chrome	16 - 22"	Right	1017237-NLA
N/A	T944A, TH94A Chrome	16 - 22"	Left	1017238-NLA
N/A	TH94A Black	16 - 22"	Right	1124235
N/A	TH94A Black	16 - 22"	Left	1124236
N/A	TH94S Black	14 - 15"	Right	1124237
N/A	TH94S Black	14 - 15"	Left	1124238
AHS4	AHS4	14 - 15"	Right	1023905-NLA
AHS4	AHS4	14 - 15"	Left	1023906-NLA
AHL4	AHL4	16 - 22"	Right	1023907
AHL4	AHL4	16 - 22"	Left	1023908
N/A	GT94HA/T94HA/T94HAB Black	16 - 22"	Right	1110315 See Note(s) B
N/A	GT94HA/T94HA/T94HAB Black	16 - 22"	Left	1110316 See Note(s) B
N/A	T94HA/T94HAB Black	16 - 22"	Right	1110317-NLA See Note(s) B
N/A	T94HA/T94HAB Black	16 - 22"	Left	1110318-NLA See Note(s) B
N/A	T94HAB Black	18 - 22"	Right	1110344 See Note(s) B
N/A	T94HAB Black	18 - 22"	Left	1110345 See Note(s) B
N/A	T94HA/T94HAB Black	16 - 22"	Right	1110348 See Note(s) B
N/A	T94HA/T94HAB Black	16 - 22"	Left	1110349 See Note(s) B
N/A	T94HS/T94HSB Black	14 - 15"	Right	1110365 See Note(s) B
N/A	T94HS/T94HSB Black	14 - 15"	Left	1110366 See Note(s) B

Pivot Slide Tube with Footplate

Legrest Pivot Slide Tube w/ Aluminum Footplate

Prior Model #	Current Model #	Chair Width	Side	Part Number
N/A	T94HS/T94HSB Chrome	14 - 15"	Right	1110367-NLA See Note(s) B
N/A	T94HS/T94HSB Chrome	14 - 15"	Left	1110368-NLA See Note(s) B

NOTE: B - T94HAB or T94HSB w/ bolt through kit.

NOTE: Footplate and pivot slide tube are attached to legrest T94HA/T94HAB using spring buttons or bolt thru hardware.

NOTE: 2" and 4" Longer Pivot Slide Tube with Footplate Assembly is not available.

Pivot Slide Tube with Footplate

Bolt Through and Spring Button Kits

Style	Part Number	Quantity per Package	Quantity Per Footrest (pr)	Quantity Per Legrest (pr)
Spring Button Hardware Kit	1188335 See Note(s) A	2	4	4
Bolt Through Hardware Kit	1101776 See Note(s) A	2	2	2
Tool Kit	1058407 See Note(s) B	1	N/A	N/A

NOTE: A - Kit includes 2 ea. Double Spring Buttons.
B - Kit includes allen wrench, key, and open end wrench. The Bolt Thru Hardware uses the Tool Kit for installation.

Footrest Pivot Slide Tube w/ Composite Footplate

Prior Model #	Current Model #	Chair Width	Side	Part Number
H903 Chrome	T943C, TH93C, T963C Chrome	16 - 22"	Right	1017224-NLA
H903 Chrome	T943C, TH93C, T963C Chrome	16 - 22"	Left	1017225-NLA
H903 Black	T943C, TH93C, T963C, Black	16 - 22"	Right	1033059
H903 Black	T943C, TH93C, T963C, Black	16 - 22"	Left	1033060
N/A	T963C, TH93S, Black	14 - 15"	Right	1035200
N/A	T963C, TH93S, Black	14 - 15"	Left	1035201
N/A	T963C, TH93S Chrome	14 - 15"	Right	1020232-NLA
N/A	T963C, TH93S Chrome	14 - 15"	Left	1020233-NLA
N/A	T963C Chrome	16 - 22"	Right	1017224-NLA
N/A	T963C Chrome	16 - 22"	Left	1017225-NLA
N/A	T93HC/T93HCB Chrome	16 - 22"	Right	1111291-NLA See Note(s) A
N/A	T93HC/T93HCB Chrome	16 - 22"	Left	1111292-NLA See Note(s) A
N/A	T93HC/T93HCB Black	16 - 22"	Right	1111289 See Note(s) A
N/A	T93HC/T93HCB Black	16 - 22"	Left	1111290 See Note(s) A
N/A	GT93HC Black	16 - 22"	Right	1194329 See Note(s) A
N/A	GT93HC Black	16 - 22"	Left	1194330 See Note(s) A

NOTE: A - T93HC/T93HCB Pivot/Slide Tube and Footplate Assemblies are attached to Footrest Supports using Spring Buttons or Bolt Thru Hardware.

NOTE: 2" and 4" Longer Pivot Slide Tube with Footplate Assembly is not available.

Pivot Slide Tube with Footplate

Legrest Pivot Slide Tube w/ Composite Footplate

Prior Model #	Current Model #	Chair Width	Side	Part Number
H904	T944C, TH94C Chrome	16 - 22"	Right	1017233-NLA
H904	T944C, TH94C Chrome	16 - 22"	Left	1017234-NLA
AHL4	AHL4	16 - 22"	Right	1023909-NLA
AHL4	AHL4	16 - 22"	Left	1023910-NLA
N/A	GT94HC/T94HC/TH9 4CB Black	16 - 22"	Right	1110311
N/A	GT94HC/T94HC/TH9 4CB Black	16 - 22"	Left	1110312
N/A	T94HC/TH94CB Chrome	16 - 22"	Right	1110313-NLA
N/A	T94HC/TH94CB Chrome	16 - 22"	Left	1110314-NLA

Footplate & Pivot Slide Tube are attached to T94HC/T94HCB Legrest using Spring Buttons or Bolt Thru Hardware.

NOTE: 2" and 4" Longer Pivot Slide Tube with Footplate Assembly is not available.

Pivot Slide Tube w/ Tubular Footplate

Prior Model #	Current Model #	Chair Width	Side	Part Number
H903 Chrome/Black	T943S	14 - 15"	Right	1030752-NLA
H903 Chrome/Black	T943S	14 - 15"	Left	1030753-NLA
H903 Chrome/Black	T943C	16 - 22"	Right	1012263-NLA
H903 Chrome/Black	T943C	16 - 22"	Left	1012264-NLA

See Support Assembly page for Pivot/Slide Tube Attaching Hardware.

Pivot Slide Tube with Footplate

Pivot Slide Tubes w/ Extra Large Footplate

Footrests

Footrest Model Number	Chair Width	Side	Part Number
T963A, T93HA, Black	16 - 22"	Right	1031924-NLA
T963A, T93HA, Black	16 - 22"	Left	1031925-NLA
T963A, T963C, H93M, T93HA Chrome	16 - 22"	Right	1031926-NLA See Note(s) A
T963A, T963C, H93M, T93HA Chrome	16 - 22"	Left	1031927-NLA See Note(s) A
T963S	16 - 22"	Right	1031928-NLA
T963S	16 - 22"	Left	1031929-NLA
T963S	14 - 15"	Right	1031930-NLA See Note(s) A
T963S	14 - 15"	Left	1031931-NLA See Note(s) A
GT93HA/T93HA, T93HAB	16 - 22"	Right	1110344
GT93HA/T93HA, T93HAB	16 - 22"	Left	1110345
T93HA	14 - 15"	Right	1110346
T93HA	14 - 15"	Left	1110347

NOTE: A - All the Footrest model numbers this is used on may not be listed here.

NOTE: 2" and 4" Longer Pivot Slide Tube with Footplate Assembly is not available.

Pivot Slide Tube with Footplate

Pivot Slide Tube w/ Extra Large Footplate

Legrest

Legrest Model Number	Chair Width	Side	Part Number
GT94HA/T94HA Black	16 - 22"	Right	1110348 See Note(s) B
GT94HA/T94HA Black	16 - 22"	Left	1110349 See Note(s) B
T94HA-Chrome	16 - 22"	Right	1110350-NLA See Note(s) A
T94HA-Chrome	16 - 22"	Left	1110351-NLA See Note(s) A
T94HA-Black	16 - 22"	Right	1124239 See Note(s) B
T94HA-Black	16 - 22"	Left	1124240 See Note(s) B
AHS4-Chrome	14 - 15"	Right	1035454-NLA See Note(s) A
AHS4-Chrome	14 - 15"	Left	1035455-NLA See Note(s) A
AHL4-Chrome	16 - 22"	Right	1035456 See Note(s) A
AHL4-Chrome	16 - 22"	Left	1035457 See Note(s) A
SH904MS, H904, T944A/C-Chrome	14 - 15"	Right	1031932-NLA See Note(s) A
SH904MS, H904, T944A/C-Chrome	14 - 15"	Left	1031933-NLA See Note(s) A
SH904M, H904, T944A/C, TH94A-Chrome	16 - 22"	Right	1031934-NLA See Note(s) A
SH904M, H904, T944A/C, TH94A-Chrome	16 - 22"	Left	1031935-NLA See Note(s) A
NOTE: A - Chrome B - Black			
NOTE: 2" and 4" Longer Pivot Slide Tube with Footplate Assembly is not available.			