

Accessories Ocean ErgoSP/(SPXL)



SPC_062A

Accessories Ocean ErgoSP/(SPXL)

| Item Number | Note | Part Number | Description | Quantity Per | |
|-------------|------|-------------|--|--------------|----------|
| | | | | Package | Assembly |
| 1 | | AP1603155 | Splash guard - Ocean Ergo Series | 1 | 1 |
| 2 | A | AP1603142 | Standard soft seat - Ocean Ergo Series | 1 | 1 |
| 3 | | AP1603143 | Soft seat insert - Ocean Ergo Series | 1 | 1 |
| 4 | B | AP1603145 | Small soft seat - Ocean Ergo Series | 1 | 1 |
| 5 | | AP1603152 | Small soft seat insert - Ocean Ergo Series | 1 | 1 |
| 6 | | 1470081 | Chest Belt Padded w/ Stainless Hardware | 1 | 1 |
| 6 | | 1470082 | Pelvis Belt Padded w/Stainless Hardware | 1 | 1 |
| 7 | | 1470077 | Armrest lock (pair) | 1 | 1 |
| 8 | | 1535077 | Back Lateral (EA) | 1 | 1 |
| 9 | | AP1643543 | Legrest Cpl., right - Ocean / Ocean Ergo Series | 1 | 1 |
| 9 | | AP1643572 | Legrest Cpl., left - Ocean / Ocean Ergo Series | 1 | 1 |
| 10 | | AP1643573 | Thigh support, right - Ocean / Ocean Ergo Series | 1 | 1 |
| 10 | | AP1643574 | Thigh support, left - Ocean / Ocean Ergo Series | 1 | 1 |
| 11 | | AP1604136 | Comfort soft backrest - Ocean (24") Ergo | 1 | 1 |
| 11 | | AP1604016 | Comfort soft backrest - Ocean (24") Ergo XL | 1 | 1 |
| 12 | | 1558769 | Soft cushion w/cover blue - Ocean / Ocean Ergo Series | 1 | 1 |
| 13 | | AP1602078 | Collection Pan holder - Ocean Ergo Series | 1 | 1 |
| 14 | | 10230 | Collection Pan with Lid | 1 | 1 |
| 15 | | AP1601986 | Front anti-tipper - Ocean (24") Ergo (XL) | 1 | 1 |
| 16 | | 16379 | Self Propelling Wheels Kit complete including brakes | 1 | 1 |
| 17 | | AP1604015 | XL backrest kit - Ocean (24") Ergo | 1 | 1 |
| 18 | | AP1606985 | Backrest tube only -Ocean VIP Ergo (separate order needed for armrest, backrest cover, headrest) | 1 | 1 |

NOTE: A - Center key shaped oval cut out is used for standard adult size individual. 11-13/16" is the approximate widest width of the oval cut out.
 B - Center key shaped oval cut out is used for smaller and young adults. 9-13/16" is the approximate widest width of the oval cut out.