



Shrouds, Rear Top and Lower Rear, Nylon



SHEP_323A



Shrouds, Rear Top and Lower Rear, Nylon

Item Number	Note	Part Number	Description	Quantity Per			
				Package	Assembly		
1	A	See Chart	Item #1 Rear Top Shroud w/ Hardware Kit	1	1		
6	B	60124360	Kit, Shroud, Rear Lower w/o Lights	1	1		
7			Shroud, Rear Lower	1	1		
8			Decal, Invacare Medallion, Reflective - Large	1	1		
9			Reflector, Red	1	2		
10			Reflector, Ring	1	2		
11			Reflector, Red	1	2		
12			Spacer, Reflector	1	2		
13			Label, Warning Motor/SSD Gearbox Replacement Sheet, Insert Parts Service	1	1		
14			C	60124349	Kit, Shroud, Rear Lower w/ Lights	1	1
7					Shroud, Rear Lower	1	1
8					Decal, Invacare Medallion, Reflective - Large	1	1
9					Reflector, Red	1	2
13					Label, Warning Motor/SSD Gearbox Replacement Sheet, Insert Parts Service	1	1
15	Harness, Rear Lights, AVIVA FX	1			1		
	Sheet, Insert Parts Service	1			1		

NOTE: A - Includes items 2-5
 B - Includes items 7-13
 C - Includes items 7-9, 13 and 15.

Item #1 Rear Top Shroud w/ Hardware Kit

Description	Standard Black	Antracite Black	Mercury Blue	Bengal Red	Silver Star	Midnight Purple
Rear Top Shroud w/ Hardware Kit	60124343-NLA	60124344-NLA	60124345-NLA	60124346-NLA	60124347-NLA	60124348-NLA
	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A

NOTE: A - Includes items 2-5.
 Nylon Colors are No Longer Available. Select from Polycarbonate Colors.