

Quick Reference Chart

When ordering 1130666, you must order 1186230 also.

If unit was built before 01/01, then the entire Manifold/Valve Assembly Kit Part #1105552 is required for Non Platinum Concentrators. See category titled "Sound Box Assembly."

Model	Compressor Assembly	4-Way Valve	P.C. Board	Sieve Beds
3LX	1101715	1101141 See Note(s) G	1078183-NLA See Note(s) A	1101707 See Note(s) O, B
5LX Before 04/02	1101715	1101141 See Note(s) G	1078183-NLA See Note(s) A	1101707 See Note(s) O, B
5LX Platinum	1106690 See Note(s) C	1101141 See Note(s) G	1078183-NLA See Note(s) A, F	1110783-NLA See Note(s) B
5LX Platinum	1106690E-NLA See Note(s) C	1101141 See Note(s) G	1078183-NLA See Note(s) A, F	1110783-NLA See Note(s) B
5LX Platinum	1157138 See Note(s) D	1101141 See Note(s) G	1108720-NLA See Note(s) E	1110783-NLA See Note(s) B
5LX Platinum	1157138 See Note(s) D	1101141 See Note(s) G	1108720E-NLA See Note(s) E	1110783-NLA See Note(s) B
5LX02 Before 7/98	1101715	1101141 See Note(s) G	1078181-NLA See Note(s) A	1101707 See Note(s) O, B
5LX02 7/98 - 03/02	1101715	1101141 See Note(s) G	1113693 See Note(s) A	1101707 See Note(s) O, B
5LX02 Platinum	1106690 See Note(s) C	1101141 See Note(s) G	1113693 See Note(s) F, A	1161699 See Note(s) O, B
5LX02 Platinum	1157138 See Note(s) D	1101141 See Note(s) G	1108731-NLA See Note(s) A	1161699 See Note(s) O, B
5LX02 Platinum	1157138 See Note(s) D	1101141 See Note(s) G	1108731E-NLA See Note(s) A	1161699 See Note(s) O, B
5SLX02T	1101715	1101141 See Note(s) G	1079444-NLA See Note(s) A	1101707 See Note(s) B, O
6LX02 Before 7/98	1075071-NLA	1101141 See Note(s) G	1078181-NLA See Note(s) A	1101707 See Note(s) B, O
6LX02 After 6/98	1075071-NLA	1101141 See Note(s) G	1113693 See Note(s) A	1101707 See Note(s) B, O
6LX02T	1075071-NLA	1101141 See Note(s) G	1079444-NLA See Note(s) A	1101707 See Note(s) O, B
6LX02H	1075071-NLA	1101141 See Note(s) G	1078181-NLA See Note(s) A	1101707 See Note(s) O, B
10LX	1130666	1101141 See Note(s) G	1128737-NLA	1144034 See Note(s) B
10LX	1130666E See Note(s) L	1101141 See Note(s) G	1128737-NLA	1144034 See Note(s) B
10LX02	1130666	1101141 See Note(s) G	1128742-NLA See Note(s) M	1144034 See Note(s) B
10LX02	1130666	1101141	1194797 See Note(s) M	1144034
10LX02	1130666E See Note(s) L	1101141 See Note(s) G	1128742-NLA See Note(s) M	1144034 See Note(s) B
Perfecto2 (IRC5P), PerfectoW (IRCPW) w/o SensO2	1162038 See Note(s) H	1101141	1143475	1161257 See Note(s) I

Quick Reference Chart

When ordering 1130666, you must order 1186230 also.

If unit was built before 01/01, then the entire Manifold/Valve Assembly Kit Part #1105552 is required for Non Platinum Concentrators. See category titled "Sound Box Assembly."

Model	Compressor Assembly	4-Way Valve	P.C. Board	Sieve Beds
Perfecto2 (IRC5PO2), PerfectoW (IRCPO2W) w/ SensO2	1162038 See Note(s) H	1101141	1143486	1161257 See Note(s) I
Perfecto2 (IRC5PO2), PerfectoW (IRCPO2W) w/ SensO2	1162038 See Note(s) H	1101141	1143486E	1161257 See Note(s) I
Perfecto2V (IRC5PO2V) w/ SensO2 See Note(s) O	1165097 See Note(s) J	1101141	1164808-NLA See Note(s) N	1165099 See Note(s) K
Perfecto2V (IRC5PO2V) w/ SensO2	1165097 See Note(s) J	1101141	1194768 See Note(s) N	1165099 See Note(s) K, O

NOTE: A - Includes Instruction Sheet.
 B - Includes (2) Sieve Beds, (2) Clamps, (2) Hose Clamps and (4) Tie Wraps.
 C - Before 4/14/04.
 D - After 4/13/04.
 E - After 11/02. Part #1108720 is NLA. See part #1108720E as a replacement.
 F - Before 12/02.
 G - If unit was built before 01/01, then the entire Manifold/Valve Assembly Kit Part #1105552 is required for Non Platinum Concentrators. See category titled "Sound Box Assembly."
 H - When replacing the compressor assembly, part number 1131250 Capacitor w/ Tie Wraps is required also.
 I - Includes (2) Sieve Beds, (2) Tie Wraps 8" Long, (4) Tie Wraps 17-1/2" Long
 J - When replacing the compressor assembly, part number 1165098 Capacitor w/ Tie Wraps is required also.
 K - Includes (2) Sieve Beds, (4) Tie Wraps 8" Long, (4) Tie Wraps 17-1/2" Long
 L - When ordering 1130666E, you must order Kit part number 1186230.
 M - 1128742 was obsoleted 5/24/17 and replaced with 1194797 (production effectivity 11/18/16). Service stock of 1128742 can only be used for units built before 11/18/16. The new 1194797 can be used for all.
 N - 1164808 was replaced with 1194768 (production effectivity 10/13/16). Service stock of 1164808 can only be used for units built before 10/13/16. New 1194768 can be used for all going forward.
 O - May contain reactivated sieve.

Quick Reference Chart

When ordering 1130666, you must order 1186230 also.

Common Replacement Filters

Part Number	Description
P0750X030	Bacteria Filter (All Models)
2000489	Cabinet Filter (All Models except Platinum 5)
1107412	Cabinet Filter (Platinum 5 Models and 10 Liter Model)
1143492	Cabinet Filter (Perfecto2/Perfecto2V/PerfectoW)
1131249	Inlet Filter (All Models)