


# Fender w/ Hardware, Drive Wheel, Nylon



## Fender w/ Hardware, Drive Wheel, Nylon

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	Item #1 Fender w/ Hardware, Drive Wheel Kit (Pair)	2	1
6	A	60124351	Kit, Fender Brackets Right and Left (Pair)	2	1
NOTE: A - NOTE: The top two Motor Mounting Screws that attach the Fender Bracket are longer than the two bottom Motor Mounting Screws.					

### Item #1 Fender w/ Hardware, Drive Wheel Kit (Pair)

Description	Standard Black	Anthracite Black	Mercury Blue	Bengal Red	Silver Star	Midnight Purple
Fender w/ Hardware, Drive Wheel Kit (Pair)	60124352-NLA See Note(s) A	60124353-NLA See Note(s) A	60124354-NLA See Note(s) A	60124355-NLA See Note(s) A	60124356-NLA See Note(s) A	60124357-NLA See Note(s) A
NOTE: A - Includes items 2-5 to accommodate both the Right and Left side. Nylon Colors are No Longer Available. Select from Polycarbonate Colors.						