


# Center Mount Rigging (CTMT-FJ/CTMT-F)



## Center Mount Rigging (CTMT-FJ/CTMT-F)

| Item Number | Note | Part Number | Description                                                                                          | Quantity Per |          |
|-------------|------|-------------|------------------------------------------------------------------------------------------------------|--------------|----------|
|             |      |             |                                                                                                      | Package      | Assembly |
| 1           | A    | 60105180    | Kit, Center Mount Front Rigging, Junior (CTMTJ/CTMT-FJ) ASBA/ADJASBA (9-1/2" - 11" Range)            | 1            | 1        |
| 2           |      |             | Rack, Center Mount Rigging Angle Adjustable                                                          | 1            | 1        |
| 3           |      |             | Support, Footrest, Center Mount, Black - Short                                                       | 1            | 1        |
| 4           |      |             | Bearing, Flange (1/2 x 13/16 x 1/2")                                                                 | 1            | 2        |
| 5           |      |             | Cleat, Button Head                                                                                   | 1            | 2        |
| 6           |      |             | Screw, Socket, Button Head w/ Patch (5/16-18 x 3/4")                                                 | 1            | 6        |
| 7           |      |             | Tube, Extension, Footrest, Short, Black - Right                                                      | 1            | 1        |
| 7           |      |             | Tube, Extension, Footrest, Short, Black - Left                                                       | 1            | 1        |
| 8           |      |             | Screw, Shoulder Hex Socket Head (5/16" x 2-1/2", 1/4-20)                                             | 1            | 1        |
| 9           |      |             | Washer (1/4 x 1/2 x 1/16")                                                                           | 1            | 1        |
| 10          |      |             | Locknut (1/4-20)                                                                                     | 1            | 1        |
| 11          |      |             | NLA - Washer, Wave (17/32 x 3/4 x 1/64")                                                             | 1            | 2        |
| 12          |      |             | Screw, Socket Button Head (1/4-20 x 1-1/2")                                                          | 1            | 2        |
| 13          |      |             | Washer, Nylon (5/16 x 11/16 x 1/16")                                                                 | 1            | 2        |
| 14          |      |             | Washer (1/2 x 7/8 x 3/64")                                                                           | 1            | 2        |
| 15          |      |             | Spacer, Nylon (1/2 x 3/4 x 1/4")                                                                     | 1            | 8        |
| 1           | A    | 60105182    | Kit, Center Mount Front Rigging (CTMT) (CTMT-F), ASBA/ADJASBA (11" - 15" Range)                      | 1            | 1        |
| 2           |      |             | Rack, Center Mount Rigging Angle Adjustable                                                          | 1            | 1        |
| 3           |      |             | Support, Footrest, Center Mount, Black - Regular                                                     | 1            | 1        |
| 4           |      |             | Bearing, Flange (1/2 x 13/16 x 1/2")                                                                 | 1            | 2        |
| 5           |      |             | Cleat, Button Head                                                                                   | 1            | 2        |
| 6           |      |             | Screw, Socket, Button Head w/ Patch (5/16-18 x 3/4")                                                 | 1            | 6        |
| 7           |      |             | Tube, Extension, Footrest, Medium, Black - Right                                                     | 1            | 1        |
| 7           |      |             | Tube, Extension, Footrest, Medium, Black - Left                                                      | 1            | 1        |
| 8           |      |             | Screw, Shoulder Hex Socket Head (5/16" x 2-1/2", 1/4-20)                                             | 1            | 1        |
| 9           |      |             | Washer (1/4 x 1/2 x 1/16")                                                                           | 1            | 1        |
| 10          |      |             | Locknut (1/4-20)                                                                                     | 1            | 1        |
| 11          |      |             | NLA - Washer, Wave (17/32 x 3/4 x 1/64")                                                             | 1            | 2        |
| 12          |      |             | Screw, Socket Button Head (1/4-20 x 1-1/2")                                                          | 1            | 2        |
| 13          |      |             | Washer, Nylon (5/16 x 11/16 x 1/16")                                                                 | 1            | 2        |
| 14          |      |             | Washer (1/2 x 7/8 x 3/64")                                                                           | 1            | 2        |
| 15          |      |             | Spacer, Nylon (1/2 x 3/4 x 1/4")                                                                     | 1            | 8        |
| 16          | B    | 60105181    | Kit, Center Mount Front Rigging w/Calf Pad - Short (CTMTJ/CTMT-FJ) ASBA/ADJASBA (9-1/2" - 11" Range) | 1            | 1        |
| 2           |      |             | Rack, Center Mount Rigging Angle Adjustable                                                          | 1            | 1        |
| 3           |      |             | Support, Footrest, Center Mount, Black - Short                                                       | 1            | 1        |
| 4           |      |             | Bearing, Flange (1/2 x 13/16 x 1/2")                                                                 | 1            | 2        |
| 5           |      |             | Cleat, Button Head                                                                                   | 1            | 2        |
| 6           |      |             | Screw, Socket, Button Head w/ Patch (5/16-18 x 3/4")                                                 | 1            | 6        |
| 7           |      |             | Tube, Extension, Footrest, Short, Black - Right                                                      | 1            | 1        |

## Center Mount Rigging (CTMT-FJ/CTMT-F)

| Item Number | Note | Part Number | Description                                                                                         | Quantity Per |          |
|-------------|------|-------------|-----------------------------------------------------------------------------------------------------|--------------|----------|
|             |      |             |                                                                                                     | Package      | Assembly |
| 7           |      |             | Tube, Extension, Footrest, Short, Black - Left                                                      | 1            | 1        |
| 8           |      |             | Screw, Shoulder Hex Socket Head (5/16" x 3", 1/4-20)                                                | 1            | 1        |
| 9           |      |             | Washer (1/4 x 1/2 x 1/16")                                                                          | 1            | 1        |
| 10          |      |             | Locknut (1/4-20)                                                                                    | 1            | 1        |
| 11          |      |             | NLA - Washer, Wave (17/32 x 3/4 x 1/64")                                                            | 1            | 2        |
| 12          |      |             | Screw, Socket Button Head (1/4-20 x 1-1/2")                                                         | 1            | 2        |
| 13          |      |             | Washer, Nylon (5/16 x 11/16 x 1/16")                                                                | 1            | 2        |
| 14          |      |             | Washer (1/2 x 7/8 x 3/64")                                                                          | 1            | 2        |
| 15          |      |             | Spacer, Nylon (1/2 x 3/4 x 1/4")                                                                    | 1            | 8        |
| 17          |      |             | Screw, Socket Button Head w/ Patch (1/4-20 x 3/4")                                                  | 1            | 6        |
| 18          |      |             | Washer, Nylon (5/16 x 7/8 x 1/8")                                                                   | 1            | 2        |
| 19          |      |             | Bracket, Calf Pad, Black                                                                            | 1            | 1        |
| 20          |      |             | Nut, T (5/16-18)                                                                                    | 1            | 2        |
| 21          |      |             | Adjuster, Calf Pad                                                                                  | 1            | 2        |
| 22          |      |             | Bracket, Mounting, Calf Pad                                                                         | 1            | 2        |
| 23          |      |             | Locknut (1/4-20)                                                                                    | 1            | 2        |
| 24          |      |             | Screw, Socket, Button Head w/ Patch (5/16-18 x 1/2")                                                | 1            | 2        |
| 25          |      |             | Pad, Calf (4-1/2" W x 7-1/4" L)                                                                     | 1            | 2        |
| 16          | B    | 60105183    | Kit, Center Mount Front Rigging w/Calf Pad - Medium (CTMT) (CTMT-F), ASBA/ADJASBA (11" - 15" Range) | 1            | 1        |
| 2           |      |             | Rack, Center Mount Rigging Angle Adjustable                                                         | 1            | 1        |
| 3           |      |             | Support, Footrest, Center Mount, Black - Regular                                                    | 1            | 1        |
| 4           |      |             | Bearing, Flange (1/2 x 13/16 x 1/2")                                                                | 1            | 2        |
| 5           |      |             | Cleat, Button Head                                                                                  | 1            | 2        |
| 6           |      |             | Screw, Socket, Button Head w/ Patch (5/16-18 x 3/4")                                                | 1            | 6        |
| 7           |      |             | Tube, Extension, Footrest, Medium, Black - Right                                                    | 1            | 1        |
| 7           |      |             | Tube, Extension, Footrest, Medium, Black - Left                                                     | 1            | 1        |
| 8           |      |             | Screw, Shoulder Hex Socket Head (5/16" x 3", 1/4-20)                                                | 1            | 1        |
| 9           |      |             | Washer (1/4 x 1/2 x 1/16")                                                                          | 1            | 1        |
| 10          |      |             | Locknut (1/4-20)                                                                                    | 1            | 1        |
| 11          |      |             | NLA - Washer, Wave (17/32 x 3/4 x 1/64")                                                            | 1            | 2        |
| 12          |      |             | Screw, Socket Button Head (1/4-20 x 1-1/2")                                                         | 1            | 2        |
| 13          |      |             | Washer, Nylon (5/16 x 11/16 x 1/16")                                                                | 1            | 2        |
| 14          |      |             | Washer (1/2 x 7/8 x 3/64")                                                                          | 1            | 2        |
| 15          |      |             | Spacer, Nylon (1/2 x 3/4 x 1/4")                                                                    | 1            | 8        |
| 17          |      |             | Screw, Socket Button Head w/ Patch (1/4-20 x 3/4")                                                  | 1            | 6        |
| 18          |      |             | Washer, Nylon (5/16 x 7/8 x 1/8")                                                                   | 1            | 2        |
| 19          |      |             | Bracket, Calf Pad, Black                                                                            | 1            | 1        |
| 20          |      |             | Nut, T (5/16-18)                                                                                    | 1            | 2        |
| 21          |      |             | Adjuster, Calf Pad                                                                                  | 1            | 2        |
| 22          |      |             | Bracket, Mounting, Calf Pad                                                                         | 1            | 2        |
| 23          |      |             | Locknut (1/4-20)                                                                                    | 1            | 2        |

## Center Mount Rigging (CTMT-FJ/CTMT-F)

| Item Number                                                                                                                                                                                                                               | Note | Part Number | Description                                          | Quantity Per |          |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------|-------------|------------------------------------------------------|--------------|----------|
|                                                                                                                                                                                                                                           |      |             |                                                      | Package      | Assembly |
| 24                                                                                                                                                                                                                                        |      |             | Screw, Socket, Button Head w/ Patch (5/16-18 x 1/2") | 1            | 2        |
| 25                                                                                                                                                                                                                                        |      |             | Pad, Calf (4-1/2" W x 7-1/4" L)                      | 1            | 2        |
| 26                                                                                                                                                                                                                                        |      | See Chart   | #26 Extension Tube Kits                              | 1            | 1        |
| 27                                                                                                                                                                                                                                        |      | See Chart   | #27 Footplate Kits                                   | 1            | 1        |
| 30                                                                                                                                                                                                                                        | C    | 60102816    | Kit, Calf Pad and Adjusting Hardware (Pair)          | 1            | 1        |
| 17                                                                                                                                                                                                                                        |      |             | Screw, Socket Button Head w/ Patch (1/4-20 x 3/4")   | 1            | 2        |
| 20                                                                                                                                                                                                                                        |      |             | Nut, T (5/16-18)                                     | 1            | 2        |
| 21                                                                                                                                                                                                                                        |      |             | Adjuster, Calf Pad                                   | 1            | 2        |
| 22                                                                                                                                                                                                                                        |      |             | Bracket, Mounting, Calf Pad                          | 1            | 2        |
| 23                                                                                                                                                                                                                                        |      |             | Locknut (1/4-20)                                     | 1            | 2        |
| 24                                                                                                                                                                                                                                        |      |             | Screw, Socket, Button Head w/ Patch (5/16-18 x 1/2") | 1            | 2        |
| 25                                                                                                                                                                                                                                        |      |             | Pad, Calf (4-1/2" W x 7-1/4" L)                      | 1            | 2        |
| 31                                                                                                                                                                                                                                        | D    | 1163374     | Kit, Footrest Pad - Small (Pair)                     | 1            | 1        |
| 29                                                                                                                                                                                                                                        |      |             | Pad, Footplate, Small - Right                        | 1            | 1        |
| 29                                                                                                                                                                                                                                        |      |             | Pad, Footplate, Small - Left                         | 1            | 1        |
| 31                                                                                                                                                                                                                                        | D    | 1163375     | Kit, Footrest Pad - Medium (Pair)                    | 1            | 1        |
| 29                                                                                                                                                                                                                                        |      |             | Pad, Footplate, Medium - Right                       | 1            | 1        |
| 29                                                                                                                                                                                                                                        |      |             | Pad, Footplate, Medium - Left                        | 1            | 1        |
| 31                                                                                                                                                                                                                                        | D    | 1163376     | Kit, Footrest Pad - Large (Pair)                     | 1            | 1        |
| 29                                                                                                                                                                                                                                        |      |             | Pad, Footplate, Large - Right                        | 1            | 1        |
| 29                                                                                                                                                                                                                                        |      |             | Pad, Footplate, Large - Left                         | 1            | 1        |
|                                                                                                                                                                                                                                           | E    | 1165505     | Kit, Hardware, Calf Pad Adjustment                   | 1            | 1        |
| 17                                                                                                                                                                                                                                        |      |             | Screw, Socket Button Head w/ Patch (1/4-20 x 3/4")   | 1            | 2        |
| 20                                                                                                                                                                                                                                        |      |             | Nut, T (5/16-18)                                     | 1            | 2        |
| 21                                                                                                                                                                                                                                        |      |             | Adjuster, Calf Pad                                   | 1            | 2        |
| 22                                                                                                                                                                                                                                        |      |             | Bracket, Mounting, Calf Pad                          | 1            | 2        |
| 23                                                                                                                                                                                                                                        |      |             | Locknut (1/4-20)                                     | 1            | 2        |
| 24                                                                                                                                                                                                                                        |      |             | Screw, Socket, Button Head w/ Patch (5/16-18 x 1/2") | 1            | 2        |
|                                                                                                                                                                                                                                           | F    | 60102817    | Kit, Bracket, Calf Pad Mounting w/ Hardware          | 1            | 1        |
| 8                                                                                                                                                                                                                                         |      |             | Screw, Shoulder Hex Socket Head (5/16" x 3", 1/4-20) | 1            | 1        |
| 9                                                                                                                                                                                                                                         |      |             | Washer (1/4 x 1/2 x 1/16")                           | 1            | 1        |
| 10                                                                                                                                                                                                                                        |      |             | Locknut (1/4-20)                                     | 1            | 1        |
| 18                                                                                                                                                                                                                                        |      |             | Washer, Nylon (5/16 x 7/8 x 1/8")                    | 1            | 2        |
| 19                                                                                                                                                                                                                                        |      |             | Bracket, Calf Pad, Black                             | 1            | 1        |
| NOTE: A - Includes items 2-15.<br>B - Includes items 2-15 and 17-25.<br>C - Includes items 17, 20-25.<br>D - Includes a pair of item 29.<br>E - Not shown. Includes items 17 and 20-24.<br>F - Not shown. Includes items 8-10, 18 and 19. |      |             |                                                      |              |          |

### #26 Extension Tube Kits

# Center Mount Rigging (CTMT-FJ/CTMT-F)

## #26 Extension Tube Kits

| Front Rigging Side                                                                                      | Short<br>9-1/4" - 12" Range | Medium<br>12-1/4" - 15" Range |
|---------------------------------------------------------------------------------------------------------|-----------------------------|-------------------------------|
| Right                                                                                                   | 1170677<br>See Note(s) A    | 1170679<br>See Note(s) A      |
| Left                                                                                                    | 1170678<br>See Note(s) A    | 1170680<br>See Note(s) A      |
| NOTE: A - Includes items 5-7. The Extension Tube Kits are used for adjustable width footplate hardware. |                             |                               |

## #27 Footplate Kits

| Footplate Side                                                                                                                                                     | Small<br>5-1/4" W x 7-1/2" L | Medium<br>5-1/2" W x 9-1/2" L | Large<br>5-7/8" W x 11-1/2" L |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------|-------------------------------|-------------------------------|
| Right                                                                                                                                                              | 1170683<br>See Note(s) A     | 1170685<br>See Note(s) A      | 1170687<br>See Note(s) A      |
| Left                                                                                                                                                               | 1170684<br>See Note(s) A     | 1170686<br>See Note(s) A      | 1170688<br>See Note(s) A      |
| NOTE: A - Includes items 10-18 and the Adjustable Width Footplate Instruction Sheet (Not Shown). The footplate mounting hardware included is for adjustable width. |                              |                               |                               |