

Front Riggings Reference Chart

70NHD (Sector Style)

NOTE: Pivot/Slide Tube attaching hardware not included.

NOTE: Pivot/Slide Tube and Footplates do not include Heel Loops.

1651 = Composite footplate with 9-7/8" curved pivot/ slide tube.
AT5543 = Adjustable angle footplate with 9-7/8" straight pivot/ slide tube.
1350 = Extra large footplate with "L" style pivot/ slide tube.

Support Assemblies - 70NHD

Rigging Type	70NHD
Right	8881053149 See Note(s) A
Left	8881053148 See Note(s) A

NOTE: A - Specify 24P frame color when ordering.

Pivot/ Slide Tube and Footplates For 70NHD

Footplate Type	1651 16" W	1651 18" and Wider	1350 18" and Wider	AT5543 ALL Widths
Right	1057869	1057879	60104886	1062756
Left	1057870	1057880	60104887	1062757
Right 2" Longer	1057871	1057881	60104888	1062760
Left 2" Longer	1057872	1057882	60104889	1062761
Right 4" Longer	1058421	1058423	60108341	1121140
Left 4" Longer	1058422	1058424	60108342	1121141