


# G-Trac-LX - LiNX Gyro Module with Mounting Hardware

Partial view of AVIVA FX Base Assembly,  
Power Module I20AL and Battery  
Harness shown for reference only



SHEP\_338A

# G-Trac-LX - LiNX Gyro Module with Mounting Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124376	Kit, Module, Gyro w/ Hardware	1	1
2			Module, GYRO LiNX (310mm)	1	1
3			Nut (M6 x 1mm)	1	1

NOTE: A - Includes items 2 and 3