


# Crutch and Cane Holder - CSCC



## Crutch and Cane Holder - CSCC

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60108349	Kit, Crutch and Cane Holder Assembly, Black - Captain Seat	1	1
2			Bracket, Crutch and Cane Holder, Black	1	1
3			Cup, Crutch/Cane Holder, Black	1	1
4			Crutch and Cane Holder Strap Assembly	1	1
5			Plug Button, Black (7/8")	1	1
6			Screw w/ Washer (#10 x 7/8")	1	2

NOTE: A - Includes items 2-6