


C629 Junction Box, Shroud Cover and Hardware



C629 Junction Box, Shroud Cover and Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124982	Kit, Junction Box, C629	1	1
2	B	60102487	Kit, Seat Frame Shroud Assembly	1	1
3	C	60124849	Kit, C629 Junction Box Extension Cable, Long 1010mm 39-49/64"	1	1
			C629 Junction Box Extension Cable, Long 1010mm 39-49/64"	1	1
			Sheet, Insert Parts Service	1	1
4	D	60124879	Kit, C629 Junction Box Extension Cable 400mm 15-3/4"	1	1
			C629 Junction Box Extension Cable 400mm 15-3/4"	1	1
			Sheet, Insert Parts Service	1	1
NOTE: A - C629 Junction Box includes the Hook and Loop Adhesive Strips. B - Use existing hardware when attaching the Seat Frame Shroud. C - Long Extension Cable used with Elevate Option. D - Extension Cable used with Non-Elevate Option.					