


# Headrest Extension Kit, Reclining Rehab Back



SHEP\_357A

## Headrest Extension Kit, Reclining Rehab Back

| Item Number                  | Note | Part Number | Description                                   | Quantity Per |          |
|------------------------------|------|-------------|-----------------------------------------------|--------------|----------|
|                              |      |             |                                               | Package      | Assembly |
| 1                            | A    | 60125832    | Kit, Headrest Extension, Reclining Rehab Back | 1            | 1        |
| 2                            |      |             | Shim                                          | 1            | 1        |
| 3                            |      |             | Plate, Headrest Extension                     | 1            | 1        |
| 4                            |      |             | Screw, Button Socket Cap (M6-1.0 x 20mm)      | 1            | 2        |
| 5                            |      |             | Locknut (M6)                                  | 1            | 2        |
| NOTE: A - Includes items 2-5 |      |             |                                               |              |          |