


LiNX CONX4 4-way Expansion Block, Shroud Cover and Hardware



LiNX CONX4 4-way Expansion Block, Shroud Cover and Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102489	Kit, LiNX CONX4 and Mounting Hardware	1	1
2		60102473	Kit, LiNX CONX4 4-way Expansion Block	1	1
3	B,A	60102487	Kit, Seat Frame Shroud Assembly	1	1
<p>NOTE: A - Use existing hardware when attaching the Seat Frame Shroud. A Foam Piece was added to the right and left side of the Seat Frame Shroud Assembly and is backward compatible to all models manufactured before the change.</p> <p>B - A Foam Piece was added to the right and left side of the Seat Frame Shroud Assembly. The addition of the two foam pieces attached to the Shroud is backward compatible.</p>					