


Drive Unit



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1592962	Handle Complete, Smoov O10	1	--
2		1592964	Tire, Smoov O10	1	--
3		1592973	Torque Receiver, Smoov O10	1	--
4		1521443	Locking Nut (M10 X 1)	1	--
5		1565715	Covering Hub Axle Nut (Z10 / Z20 / P15)	1	--
	A	1456190	Screw, Cylinder Head (DIN 912 - M 6 x 16 - 8.8)	1	--
	A	1592841	Cover, Charging Socket O10	1	--

NOTE: A - Not Shown