


Seat Upholstery



Seat Upholstery
Sail Cloth

Seat Upholstery - Sail Cloth - Non Heavy Duty

Chair Width	14" Deep	15" Deep	16" Deep	17" Deep	18" Deep	19" Deep	20" Deep
14"	1032080	1032081	1032082	1032083	1032084	N/A	N/A
15"	1032085	1032086	1032087	1032088	1032089	N/A	N/A
16"	1032090	1032091	1032092	1032093	1032094	N/A	N/A
17"	1032095	1032096	1032097	1032098	1032099	N/A	N/A
18"	1032100	1032101	1032102	1032103	1032104	1111857	1111858
19"	N/A	N/A	1032105	1032106	1032107	1111859	1111860
20"	N/A	N/A	1032108	1032109	1032110	1111861	1111862

Seat Upholstery

Seat Upholstery - Heavy Duty

Chair Width	18" Deep	19" Deep	20" Deep	21" Deep	22" Deep
18"	1107534	1111869	1111838	N/A	N/A
19"	1107537	1111839	1111840	N/A	N/A
20"	1107540	1107541	1107542	1107543	1107544-NLA
22"	60125703	60125704	60125705	60125706	1107545

Seat Upholstery Hardware Kit

Chair Style	Hardware Kit	Screw Quantity/Size
Non Heavy Duty	1046294-NLA	16 ea. (#10-32 x 3/4")
Heavy Duty	1076845	18 ea. (#10-16 x 7/8")