


Lever Assembly for Leg Spread ISA COMPACT



Lever Assembly for Leg Spread ISA COMPACT

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	AP1651334-701 6	Lever Assembly for Manual Leg Spread	1	1
2		SP1654296	Lever, Manual Leg Spread	1	1
3		SP1645498	Plate, Leg Spread Lever	1	1
4		SP1449580	Plug, Plastic, Leg Spread Lever	1	1

NOTE: A - Includes items 2-4