




# Compressor Assembly, Platinum 5NXG



SRE\_042A



Assembled view shown for reference only

## Compressor Assembly, Platinum 5NXG

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60127170	Kit, Compressor Assembly, Platinum 5NXG 120V AC, NA	1	1
2			Compressor - Platinum 5NXG 120V AC, NA	1	1
3			Spring, Compression, Platinum 5NXG	1	4
4			Screw, Socket Head Cap w/ Patch (3/8-16 x 5/8")	1	4
5			Cup, Compressor Spring, Platinum 5NXG	1	4
6			Tube, Molded Compressor	1	1
7			Clamp, Hose	1	1
8			Foam, Compressor	1	2
9			Insulator, Compressor Spring	1	1
			Sheet, Insert Parts Service	1	1
10	B	60127172	Kit, Tube, Heat Exchanger to Compressor	1	1
6			Tube, Molded Compressor	1	1
7			Clamp, Hose	1	1
			Sheet, Insert Parts Service	1	1

NOTE: A - Includes items 2-9. For North America (NA) models  
 B - Includes items 6 and 7.