


Inlet Filter Assembly



Platinum 5NXG Lower Assembly
shown for reference only

Inlet Filter Assembly

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60127147	Kit, Inlet Sound Filter, Platinum 5NXG	1	1
			Filter, Inlet Sound Reduced	1	1
			Sheet, Insert Parts Service	1	1
2	A	60127148	Kit, Inlet Sound Filter Flange and Tubing, Platinum 5NXG	1	1
3			Flange, Inlet Filter, Platinum 5NXG	1	1
4			Tube, Inlet to Compressor Molded, Platinum 5NXG	1	1
5			Washer, Inlet Filter	1	1
			Sheet, Insert Parts Service	1	1
5		60127175	Kit, Inlet Filter Shell Cover	1	1
			Cover, Inlet Filter Shell	1	1
			Sheet, Insert Parts Service	1	1

NOTE: A - Includes items 3-5.