


Fans, Compressor



Platinum 5NXG Base Assembly
and Shell Component shown
for reference only

NA Power Cord Assembly
shown for reference only

SRE_048A

Fans, Compressor

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60127169	Kit, Dual Compressor Fans, 12VDC	1	1
2			Fan, Compressor, 12VDC	1	2
3			Guard, Fan	1	2
4			Grommet, Fan Guard	1	8
			Sheet, Insert Parts Service	1	1

NOTE: A - Includes items 2-4