




Accessories for Model 7740, 7740P, 7740A (Octagon Tube)




**7950 and 7950P
Upper Bed Bracket**




**7953 and 7953P
Lower Bed Bracket
with Wall Bumper**




**7954 and 7954P
Lower Bed Bracket**



**07800X001 and 1125219
Hinge Clamp**



**7729 and 7729P
Swivel Clamp
and Bracket**



**7941 and 7941P
Grab Bar, Chain
and Clamp**

Accessories for Model 7740, 7740P, 7740A (Octagon Tube)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		10301X000-NLA	NLA - Hook, "S"	1	2
2	H,B	17571X000-NLA	NLA - Bracket, Upper, Chrome (Includes Clamp and Wing Nut)	1	1
2	C	7729P-NLA	NLA - Bracket, Upper - Grey (Includes Clamp and Wing Nut)	1	1
3		17572X000	Rod, Clamp	1	1
4		16324X000	Nut, Wing (3/8"-16)	1	2
5		14130X000-NLA	NLA - Plunger	1	1
6	B	1009042-NLA	NLA - Clamp, Swivel w/ Plunger and Wing Knob, Chrome	1	1
6	C	1137116-NLA	NLA - Clamp, Swivel with Plunger and Wing Knob, Grey	1	1
7		17129X000-NLA	NLA - Knob, Wing	1	1
8	E,B	17074X000-NLA	NLA - Bracket, Trapeze, Chrome	1	1
8	C	1135738-NLA	NLA - Bracket, Trapeze, Grey	1	1
9	A	17077X000-NLA	NLA - Hook Rod	1	2
10		40186X030	Tip, Rubber, Grey (1")	1	1
11	D,B	17061X000-NLA	NLA - Plate, Back with Clamp, Chrome	1	1
11	C	1135728-NLA	NLA - Plate, Back with Clamp, Grey	1	1
12		17127X000	Pin, Spiral (5/32" x 1-3/4")	1	1
13		17128X000-NLA	NLA - Pin, Spiral (5/16" x 2-1/2")	1	1
14		1008630-NLA	NLA - Stud, Threaded (3/8-16 x 2-1/4")	1	1
15	G,B	17121X000-NLA	NLA - Clamp, Hinge, Chrome	1	1
15	J	1135725-NLA	NLA - Clamp, Hinge, Grey	1	1
16	G,B	17115X000-NLA	NLA - Clamp, Hinge, Slotted Side, Chrome	1	1
16		1135733-NLA	NLA - Clamp, Hinge, Slotted Side, Grey	1	1
17	B	17101X000-NLA	NLA - Grab Bar, Trapeze, White	1	1
17	C	1109130-NLA	NLA - Grab Bar, Trapeze, Grey	1	1
18		95103X000-NLA	NLA - Chain, 8 Link	1	1
19	F,B	17078X000-NLA	NLA - Bracket, Bed Frame, Chrome	1	1
19	C	1135737-NLA	NLA - Bracket, Bed Frame, Grey	1	1
20		1046727	Cap, Plug, Grey	1	2
21		1011755	Foam, Adhesive	1	1

NOTE: A - As of 8/98 Hook Rod is NLA, replace with Clamp Rod - item #3.
 B - Before 9/9/05.
 C - After 9/8/05.
 D - This part is No Longer Available, please use part # 7950P found in the "Universal Fracture Frame" section of this book.
 E - This part is No Longer Available, please use part # 7953P found in the "Offset Trapeze Bar" section of this book.
 F - This part is No Longer Available, please use part # 7954P found in the "Universal Fracture Frame" section of this book.
 G - This part is No Longer Available, please use part # 1125219 found in the "Accessories for Model 7030 & 7030P (Octagon Tube)" section of this book.
 H - The chrome part is NLA please replace with the grey direct replacement part.
 J - This parts is available for both the 7740, 7740P, 7740A Trapeze.