


# Fender w/ Hardware, Drive Wheel



## Fender w/ Hardware, Drive Wheel

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	Item #1 Fender w/ Hardware, Drive Wheel Kit	1	1
6		60127913	Kit, Fender Brackets (Right and Left)	1	1
7			Bracket, Fender, Black - Right	1	1
8			Bracket, Fender, Black - Left	1	1
9			Screw, Socket Head Cap (M4 x 0.7 x 12)	1	6
10			NLA - Washer, Lock Internal Tooth (#10)	1	6
			Sheet, Insert Parts Service	1	1

NOTE: A - Includes items 7-10.

### Item #1 Fender w/ Hardware, Drive Wheel Kit

Description	Obsidian Black (170P)	Ruby Red (171P)	Peacock Blue (172P)	Pink Flamingo (173P)	Purple Passion (174P)	Denim Navy (175P)
Fender w/ Hardware, Drive Wheel Kit	60129165	60129166	60129167	60129168	60129169	60129170
	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A

NOTE: A - Includes items 2-5 to accommodate one side only.

### Item #1 Fender w/ Hardware, Drive Wheel Kit (continued)

Description	Lime Twist (176P)	Artic Blue (177P)
Fender w/ Hardware, Drive Wheel Kit	60129171	60129172
	See Note(s) A	See Note(s) A

NOTE: A - Includes items 2-5 to accommodate one side only.