



Electronics - GYRO Module and Mounting Hardware



Rear Side of AVIVA Storm RX
Base Assembly shown for
reference only



LiNX Power Module 120AL
Shown for reference only

Electronics - GYRO Module and Mounting Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60127936	Kit, Module, GYRO and Mounting Hardware	1	1
2			Module, GYRO LiNX (310mm)	1	1
3			Nut (M6 x 1mm)	1	1
			Sheet, Insert Parts Service	1	1

NOTE: A - Includes items 2 and 3