


Rear Wheel Hardware



Rear Wheel Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	47203X024-NLA	NLA - Bolt, Hex Head (7/16-20 x 4")	1	1
1	B	47203X025	Bolt, Hex Head Cap (7/16-20 x 4-3/4")	1	1
2	A	40157X011-NLA	NLA - Washer, Steel (7/16 x 15/16 x 1/16")	1	1
2	B	47183X037-NLA	NLA - Spacer (1/2 x 7/8 x 7/8")	1	1
3		7930352-NLA	NLA Locknut (7/16-20)	25	1
4	C	1027600-NLA	NLA - Cap, Black	1	1
NOTE: A - Permanent Arm Chairs and Non-Reclining B - Removable Arm Chairs and All Recliners C - DTM Wheel Only					