


# Aquatec Ocean Shower Commode Accessories



SPC\_076A - SPS1548776

## Aquatec Ocean Shower Commode Accessories

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		16341	Soft Seat w/ Key Cut	1	1
7		1544469	Variable soft seat cpl. - Ocean / VIP / Ergo Series	1	1
7		1545659	Clip, Variable Softseat	1	2
8		10167-10	Safety Support Bar	1	1
9		1470078	Safety Support Bar , curved	1	1
11		1470081	Chest Belt Padded w/ Stainless Hardware	1	1
11		1470082	Pelvis Belt Padded w/Stainless Hardware	1	1
12		1470077	Armrest lock (pair)	2	1
14		16379	Self Propelling Wheels Kit complete including brakes	2	1
15	A	AP1643574	Thigh support, left - Ocean / Ocean Ergo Series	1	1
15	B	AP1643573	Thigh support, right - Ocean / Ocean Ergo Series	1	1
17	C	AP1643572	Legrest Cpl., left - Ocean / Ocean Ergo Series	1	1
17	D	AP1643543	Legrest Cpl., right - Ocean / Ocean Ergo Series	1	1
18		1535077	Back Lateral (EA)	1	1
19		1558769	Soft cushion w/cover blue - Ocean / Ocean Ergo Series	1	1
20		10230	Collection Pan with Lid	1	1
22		1575299	Ocean Contour Hygienic Fixed Back Support	1	1
NOTE: A - Replaces 1556894 B - Replaces 1556893 C - Replaces 1546705 D - Replaces 1546439  Items not listed are non-saleable					